




29 Wells Avenue Yonkers, NY 10701 914-963-2040 800-431-2972

SHEEPSKIN WHEELCHAIR ACCESSORIES (Armrest Pads)


Reorder #	Unit
703010	Pair



LAUNDERING INSTRUCTIONS

These Pads may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.

Advanced design for nursing and therapy.™





WHEELCHAIR LEG PAD

 Re-Order #	Size	Unit
703020	12"x 24"	Each



PURPOSE

High-pile synthetic sheepskin fabric protects your patient's legs against contact with metal wheelchair parts. Effective for preventing patient's feet from slipping off footrests. Sized to fit all wheelchairs, the Leg Pad serves both as a cushion and a leg support when the footrests are in the raised position.

APPLICATION

- 1. Place pad lengthwise across the metal bars that support footrests.
- 2. With fluff side facing out, fold ends around bars.
- **3.**The pad should be positioned as shown in illustration above.
- **4.**The pad is held in place by securing the straps through the corset buckles.


LAUNDERING INSTRUCTIONS

The Wheelchair Leg Pad is made from 100 percent polyester synthetic sheepskin fabric. It can be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.





29 Wells Avenue Yonkers, NY 10701 914 963-2040 800 431-2972



WHEELCHAIR FOOTREST PADS

Re-Order #	Previous #	Unit	
703040	9035	Pair	

These super-soft synthetic sheepskin pads cushion sensitive feet against metal footrests. Designed to fit all wheelchairs, the pads are held securely in place with elastic bands.

LAUNDERING INSTRUCTIONS

The Footrest Pads may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.



Advanced Design for Nursing and Therapy.

914-963-2040 800-431-2972

WHEELCHAIR PLUSH PADS (For Seat and Backrest)

Re-Order #		Size	Unit
	703050	18"x16"	Pair (Seat and Chair)

These pads add a touch of luxury to any wheelchair. The high synthetic sheepskin provides pressure relief plus improved air circulation. Packaged as a set of two units, one for the seat and one for the backrest, each pad has tie straps that secure it to the wheelchair.

LAUNDERING INSTRUCTIONS

Ski-Care's Plush Pads may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.

Advanced design for nursing and therapy.™

