



ReOrder	Description	Size	Units
909220	Stop Strip w/Stop Sign	12" x 50"	Each
909221	Stop Strip - Yellow (no stop sign)	12" x 50"	Each
909225	Stop Strip w/Magnetic Attachment	12" x 50"	Each
909226	Stop Sign w/Velcro	12" x 12"	Each
909217	Stop Strip alarm System w/Velcro Mounts	12" x 50"	Each
909218	Stop Strip alarm System w/Velcro Mounts	18" x 50"	Each

Alarm included with model 909217

**Purpose**

These brightly colored strips help to prevent accidents by directing residents away from potential hazard areas. Promotes resident privacy. Placed across a room doorway, this brightly colored strip often stops the cognitively impaired wanderer from entering another resident's room. It also acts as a reminder for the resident inside a room not to leave without first notifying a staff member. Secured with Velcro<sup>®</sup>, the Stop-Strip is not a physical barrier, so it may be used as part of a restraint reduction program. The stop strip barrier is available with or without an alarm. Procedure A discusses how to apply the unit without the alarm whereas Procedure B covers application with the alarm.

**Easy to Use**

**Procedure A**

Peel the backing off the two Velcro<sup>®</sup> pressure-sensitive strips and attach them to the left and right sides of the door frame or to the wall on either side of the door frame. These strips should be positioned at mid height on the door frame or on the wall next to the door frame. (Fig. 1) **IMPORTANT:** Attach the pressure sensitive strips to the corridor side of the door frame or wall so that the Stop-Strip is outside the room. (Fig. 2).

**Procedure B**

1. Install 9-volt battery. Battery compartment is on side of alarm unit. The Stop Strip alarm system includes a **long** and **short** Velcro<sup>®</sup> strip. Peel the backing off the **short** strip and attach it to the **left** side of the door frame. Attach the **long** strip to the **right** side of the door frame. Strips should be attached at the resident's shoulder height (see Fig. 2).
2. Attach alarm unit to the **short** strip on the **left** side of the door (see Fig. 3).
3. Attach the Velcro<sup>®</sup> material on the door strip to the **long** strip on the **right** side of the door. Bring door strip across the door and attach the magnet sewn into the door strip to the strip on the alarm unit (see Fig. 1).
4. Set alarm volume by moving the switch to the Hi or Low position.




Fig. 1




Fig. 2




Fig. 3

**Caution**

- Test the Stop Strip alarm before each use.
- Move the switch to the Hi position and remove the magnet.
- If alarm is weak or does not sound, replace 9-volt battery.
- Do not place the Stop-Strip across fire and emergency exit doors.