

Wheelchair Guide

Model Style	SIZING					FRONT RIGGING OPTIONS		ARMREST OPTIONS					WARRANTY
	Seat Height	Back Height	Seat Width	Seat Depth	Weight Capacity	Swing-Away Footrests	Swing-Away Elevating Legrests	Fixed	Removable	Adjustable Height	Desk Length	Full Length	Warranty Period
STANDARD MODELS													
Charisma™	19"	16"	16" - 18"	16"	250 lbs	X	X	X	X		X	X	2 Year Limited
114061	19"	16"	16"	16"	250 lbs	X		X				X	2 Year Limited
114062	19"	16"	16"	16"	250 lbs		X	X				X	2 Year Limited
114063	19"	16"	16"	16"	250 lbs	X			X			X	2 Year Limited
114064	19"	16"	16"	16"	250 lbs		X		X			X	2 Year Limited
114065	19"	16"	16"	16"	250 lbs	X			X		X		2 Year Limited
114066	19"	16"	16"	16"	250 lbs		X		X		X		2 Year Limited
114067	19"	16"	18"	16"	250 lbs	X			X			X	2 Year Limited
114068	19"	16"	18"	16"	250 lbs		X		X			X	2 Year Limited
114529	17" - 19"	16"	18"	16"	250 lbs	X			X			X	2 Year Limited
114530	17" - 19"	16"	18"	16"	250 lbs		X		X			X	2 Year Limited
114531	17" - 19"	16"	18"	16"	250 lbs	X			X		X		2 Year Limited
114532	17" - 19"	16"	18"	16"	250 lbs		X		X		X		2 Year Limited
114505	17" - 19"	16"	20"	16"	250 lbs	X			X			X	2 Year Limited
114506	17" - 19"	16"	20"	16"	250 lbs		X		X			X	2 Year Limited
114507	17" - 19"	16"	20"	16"	250 lbs	X			X		X		2 Year Limited
114508	17" - 19"	16"	20"	16"	250 lbs		X		X		X		2 Year Limited
Sentra EC													
114583	17.5" - 19.5"	18"	22" - 24"	18"	450 lbs	X	X		X		X	X	3 Year Limited
114584	17.5" - 19.5"	18"	22"	18"	450 lbs	X			X		X		3 Year Limited
114585	17.5" - 19.5"	18"	22"	18"	450 lbs	X			X			X	3 Year Limited
114586	17.5" - 19.5"	18"	22"	18"	450 lbs		X		X			X	3 Year Limited
114587	17.5" - 19.5"	18"	24"	18"	450 lbs	X			X		X		3 Year Limited
114588	17.5" - 19.5"	18"	24"	18"	450 lbs		X		X		X		3 Year Limited
114589	17.5" - 19.5"	18"	24"	18"	450 lbs	X			X			X	3 Year Limited
114590	17.5" - 19.5"	18"	24"	18"	450 lbs		X		X			X	3 Year Limited
BARIATRIC MODELS													
Sentra Heavy-Duty	17.5" - 19.5"	18"	22" - 24"	18"	500 lbs	X	X		X	X	X	X	Limited Lifetime
114591	17.5" - 19.5"	18"	22"	18"	500 lbs	X			X	X	X		Limited Lifetime
114592	17.5" - 19.5"	18"	22"	18"	500 lbs		X		X	X	X		Limited Lifetime
114593	17.5" - 19.5"	18"	22"	18"	500 lbs	X			X	X		X	Limited Lifetime
114594	17.5" - 19.5"	18"	22"	18"	500 lbs		X		X	X		X	Limited Lifetime
114595	17.5" - 19.5"	18"	22"	18"	500 lbs	X			X		X		Limited Lifetime
114596	17.5" - 19.5"	18"	22"	18"	500 lbs		X		X		X		Limited Lifetime
114597	17.5" - 19.5"	18"	22"	18"	500 lbs	X			X			X	Limited Lifetime
114598	17.5" - 19.5"	18"	22"	18"	500 lbs		X		X			X	Limited Lifetime
114599	17.5" - 19.5"	18"	24"	18"	500 lbs	X			X	X	X		Limited Lifetime
114600	17.5" - 19.5"	18"	24"	18"	500 lbs		X		X	X	X		Limited Lifetime
114601	17.5" - 19.5"	18"	24"	18"	500 lbs	X			X	X		X	Limited Lifetime
114602	17.5" - 19.5"	18"	24"	18"	500 lbs		X		X	X		X	Limited Lifetime
114603	17.5" - 19.5"	18"	24"	18"	500 lbs	X			X		X		Limited Lifetime
114604	17.5" - 19.5"	18"	24"	18"	500 lbs		X		X		X		Limited Lifetime
114605	17.5" - 19.5"	18"	24"	18"	500 lbs	X			X			X	Limited Lifetime
114606	17.5" - 19.5"	18"	24"	18"	500 lbs		X		X			X	Limited Lifetime
RECLINING BACK MODELS													
Silver Sport	20.5"	21" - 31"	16", 18", 20"	16"	300 lbs	X	X		X		X	X	Limited Lifetime
114721	20.5"	21" - 31"	16"	16"	300 lbs				X		X		Limited Lifetime
114722	20.5"	21" - 31"	18"	16"	300 lbs				X			X	Limited Lifetime
114723	20.5"	21" - 31"	20"	16"	300 lbs				X		X		Limited Lifetime