



Call Crest Technical Support for assistance at: **1-800-328-8908**

www.cresthealthcare.com

©2010 Crest Healthcare Supply



Crest 7700 Repairable Momentary Call Cords

Crest's 7700 Series Call Cords use a repairable push-button pendant and can be ordered with a plug to match your nurse call system. Each call cord includes a security clip. Available in white, grey or beige.

Specifications

ELECTRICAL

Max Voltage Rating: 28VAC/DC

OXYGEN SAFETY: Meets UL 1069 Requirements

LATEX INFORMATION: This Product is Latex Free

MECHANICAL

Button Life: >1,000,000 cycles

Switch Type: Normally Open SP/ST

SIZE MEASUREMENTS (PENDANT)

Height: 2.25"

Diameter: 1" (at widest point)

WARRANTY: 1 year