



Call Crest Technical Support for assistance at: **1-800-328-8908**

www.cresthealthcare.com

©2010 Crest Healthcare Supply



Crest DuraCall™ Momentary Call Cords

DuraCall is a premium momentary call cord designed to be extremely durable, water resistant and long lasting. DuraCall is compatible with most major nurse call systems. The large, winged style pendant makes it easy and comfortable to grip. The completely sealed pushbutton has a raised button design that provides both a tactile and audible “click” to assure residents a call has been placed.

Specifications

ELECTRICAL

Max Voltage Rating: 28VAC/DC

OXYGEN SAFETY: Meets UL 1069 Requirements

LATEX INFORMATION: This Product is Latex Free

MECHANICAL

Button Life: >1,000,000 cycles

Switch Type: Normally Open SP/ST

SIZE MEASUREMENTS (PENDANT)

Height: 2.75"

Width: 1.625" (at widest point)

WARRANTY: 2 Years