

Crest HD2™ Bed Controls

The Cost-Effective and Durable Solution



4 Button Model

Controls Head and Foot functions



6 Button Model

Controls Head, Bed and Foot functions



8 Button Model

Featuring auto-contour that controls the Head and Foot simultaneously

Crest HD2 replacement bed controls are durable and comfortable to use. Low-voltage HD2 bed controls were designed to be completely sealed with universal graphics that will never flake, pull or peel off. If a cable replacement is needed, the “plug-n-play” solderless cable is simple to replace. HD2 Bed Controls are white and feature an 8-foot cord with durable strain relief.

- Lightweight and ergonomic
- Large buttons have an audible & tactile click
- Long life mechanical switches
- Sealed overlay is easy to sanitize & keep clean

Try our Bed Control Holder!

Made of durable yet flexible polypropylene which holds HD2 bed controls snugly and can be hooked on bed rails or headboards without damage concerns.

HD2 Bed Controls are available for many different bed types. Call for details.



"Crest replacement bed controls are the most durable we have found. They last four times longer than the controls we get from the OEM. Their shape is more comfortable for our residents to use, too." – Robert Johnston, Redstone Highlands, Pennsylvania

CREST
HEALTHCARE®

800.328.8908

800.369.9207

CrestHealthcare.com



Don't just take our word for it.
Here is what our customers say!

"I don't buy my bed controls from anyone else."

– Albaro, CA

"The controls are easy for the residents to see."

– James, NM

"The shape is ideal for our residents – the corners are curved, and there are no sharp edges to hurt the seniors."

– Jose, CA

"Crest has the best price on bed controls."

– Dwayne, NV

"They are built with quality and last for years. They last many times longer than other brands. If you put them side by side with others, you can just see Crest's is a better product."

– Jim, NM

"The price is right."

– Brian, PA

"Crest bed controls are very durable. They last longer than the OEM bed controls, and I don't have any issues with the boards."

– Bill, MD

"I work with ten facilities, and I tell them all to buy Crest bed controls."

– Jim, NM

"I do my own cord replacements. It's very easy and takes just a few minutes. Buying a repair kit is a lot cheaper than buying a new bed control, and replacing the cable doesn't take any extra time."

– Ray, IN



"They are easy to clean. They are sealed and easy to wipe down. There aren't buttons that stick out and trap dirt."

– Jose, CA