



Crest Healthcare

Jeron® Pro-Alert™ 640 & Provider® 680 Part List

Your guide to Jeron® Pro-Alert® 640 & Provider 680 system parts available from Crest that are durable, economical and designed with infection control in mind.

CREST
HEALTHCARE®

© 2023 Crest Healthcare



Enhancing the environment of care for those in need and those who care

Crest Healthcare is a leading provider of a wide range of healthcare products and technical services. Historically specialized in high-quality nurse call replacement parts, and today producing nurse call system components for nearly every current system, Crest has been a trusted source for innovative product development and production since 1967.

With customer satisfaction as our driving goal, Crest products provide quality, value, innovation, variety and durability. Working together to guarantee a positive experience, our in-house team includes customer service, technical support, engineering, manufacturing and quality assurance experts.

Engineering

Our engineering team is continually creating and testing new communication products. With the use of destructive and non-destructive in-house testing and prototype evaluations, Crest ensures every product provides durability and a long service life. Custom services are available to OEMs and other commercial customers.

Manufacturing

We build and assemble all the Crest products in-house and work closely with our manufacturing partners to ensure our same high quality standards are met.

Quality Assurance

We conduct testing on 100% of finished products we manufacture prior to shipment. Crest products meet or exceed demanding safety standards up to and including UL 1069 Recognition.

Quality Certifications

Crest strives to offer you high quality products to support your facility needs. Crest Healthcare is ISO 13485:2016 certified provider of Pillow Speakers, Nurse Call Stations, Nurse Call Cords, and Bed Control Equipment. Our certification represents our commitment to on-going quality improvements and a dedication to higher levels of performance, service and quality products. We received our first certification in 1999 (ISO 9001) and Crest continues to improve quality in all aspects of designing, manufacturing and customer service.



 800.328.8908
 800.369.9207
 CrestHealthcare.com

Jeron® Pro-Alert™ 640 & Provider™ 680
Nurse Call Components

Pillow Speakers

Elite Digital™4
Elite Infrared™5
Replacement Cable Assemblies6
Magnetic Break-a-Way Cable System6
Remote Call-IR™7
Slimline™8

Call Cords

Momentary Call Cords

EconoCall™9
CleanCall™ 10
DuraCall™ 11
DuraColor™ 12
DuraGlo™ 12
DuraLight™ 12
SoftCall™ 13

Specialty Call Cords

BreathCall® 14
GeriCall™ 14
PressCall® 15
PadCall® 15
TouchCall™ 16

Additional Components

Adapters..... 17
Station Savers..... 17
Call Cord Extenders..... 18
Dummy Plug 18
Receptacles..... 18

Stations

Patient Station..... 19
Pullcord Station..... 19




EliteDigital™

Crest Elite Digital™ Pillow Speakers

Designed for use with digital TVs, Autosense™ Technology automatically identifies the brand of digital healthcare TV and sends the correct signal. They can be moved to any room, regardless of healthcare television brand.

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille



Specifications	
Height:	9"
Width:	4 3⁄5"
Thickness:	1 3⁄4"
Cord Length:	8'
Case Material:	ABS, UL94V-0
Cord Connector:	Plug-and-Play (solderless)
Power Requirements:	Two AAA alkaline batteries required (included) (Infrared only)
Switch Rating:	30 VDC @ 100 mA
Warranty:	2-year warranty against defects in materials and workmanship from date of purchase



Part No.	System	Nurse Call	TV Control	Other Features	Price
XN3801XBAC60A10N1†	Jeron Pro-Alert 640 and Provider 680	Nurse	10-key		121.35
XN3801XCAC60A10N1†	Jeron Pro-Alert 640 and Provider 680	Nurse	10-key	1 light	121.35
XN3801EBAC60A10N1†	Jeron Pro-Alert 640 and Provider 680	Nurse	Up/Down		129.74
XN3801ECAC60A10N1†	Jeron Pro-Alert 640 and Provider 680	Nurse	Up/Down	1 light	129.07

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.




EliteInfrared™

Crest Elite Infrared™ Pillow Speakers

The Crest Elite Infrared Pillow Speaker is the only pillow speaker that combines nurse call and TV with external DVD or Set Top Box (STB) controls. It combines the ease of a universal remote with traditional pillow speaker functionality such as nurse call and environmental controls. It is also useful where there is no cross-room wiring or with TVs that do not have a pillow speaker port. Uses two AAA batteries (included).

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille



Specifications	
Height:	9"
Width:	4 3⁄5"
Thickness:	1 3⁄4"
Cord Length:	8'
Case Material:	ABS, UL94V-0
Cord Connector:	Plug-and-Play (solderless)
Power Requirements:	Two AAA alkaline batteries required (included) (Infrared only)
Switch Rating:	30 VDC @ 100 mA
Warranty:	2-year warranty against defects in materials and workmanship from date of purchase




Part No.	System	Nurse Call	TV Control	Other Features	Price
XT3801XBAC60Z10N1†	Jeron Pro-Alert 640 and Provider 680	Nurse	10-key		148.46
XT3801XCAC60Z10N1†	Jeron Pro-Alert 640 and Provider 680	Nurse	10-key	1 light	148.46
XT3801ECAC60Z10N1†	Jeron Pro-Alert 640 and Provider 680	Nurse	Up/Down	1 light	142.57
XT3801PCAC60Z10N1†	Jeron Pro-Alert 640 and Provider 680	Nurse	10-key	STB	172.03

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.

Replacement Cable Assemblies for Elite Pillow Speakers

Replacement pillow speaker cable assemblies for Elite Digital™ and Elite Infrared™ are easy to replace using "plug-and-play" connectors. Standard white cables are 8' with molded-in strain relief. To verify the correct cable for your pillow speaker or for replacement assistance, please contact Crest Technical Support.



Part No.	System	Color	Length	Price
XN3801W-08	Jeron Pro-Alert 640 and Provider 680	White	8'	27.95




Crest Magnetic Break-A-Way Cable System

The Magnetic Break-a-Way Cable System saves both wall stations and pillow speakers from damage due to cable stress. Any cable pressure over 6 lb will break the magnetic connection and disconnect the pillow speaker without damage to the wall station. The magnets guide the connection, ensuring the pillow speaker is correctly and quickly connected - no more fumbling or broken pins.

- For use with Crest devices
- Eliminate downtime & cost of wall station repairs
- Save on pillow speaker replacements
- Cable length is approximately 24"
- Also works with call cords

LATEX-FREE



Part No.	System	Color	Length	Price
116075	Jeron Pro-Alert 640 and Provider 680	White	2'	74.40



All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V1. Revised 07/25/23.




Crest Remote Call-IR™ TV Control & Nurse Call

Replace your call cord and remote control with one simple device. The Crest Remote Call-IR controls the TV using IR signals like a universal remote while still providing the safety of reliable wired nurse call.

Remote Call-IR plugs into your nurse call station to ensure it stays within the room and is durable enough to withstand drops and daily use. The 10-key overlay style is great for emergency rooms or private patient rooms, where the room is not wired for TV sound or the sound control of a pillow speaker is not needed. The large button overlay is perfect for kids, geriatric patients or residents with limited dexterity. Features simplified control of TV power, channel and volume with big, easy-to-see buttons. Both styles feature an on-board library of codes that can be used to program it to work with almost any TV. Uses two AA batteries (not included) that last as long as in typical TV remotes.

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille
- Bed A & B functionality
- Fits Crest Pillow Speaker Keeper



Part No.	System	Nurse Call	Style	Length	Price
XP441E0B600XXX2 †	Any station with 1/4" jack or Type B or Y adapter	Nurse	Large Button	10'	81.95
XP441X0B600XXX2 †	Any station with 1/4" jack or Type B or Y adapter	Nurse	10-key	10'	89.95

Replacement Cable for Remote Call IR

Part No.	System	Length	Price
XP4401W-10	Any station with 1/4" jack or Type B or Y adapter	10'	33.91

! Programming instructions and code list can be found at www.cresthealthcare.com under the "Resources" tab

Specifications

Height:	6 ½"
Width:	Tapers, 2 ¾" - 2 ¼"
Thickness:	1 ⅝" at thickest point
Weight:	13 ⅞ oz (with batteries)
Cord Length:	10'
Case Material:	ABS, Flame rated (UL94V-0)
Cord & Strain Relief:	Flexible molded PVC
Cord Connector:	Plug-and-Play (solderless)
Plug Style:	2-conductor or 3-conductor ¼" phone plug (call Crest for available adapters)
Power Requirements:	Two AA alkaline batteries (sold separately)
Switch Rating:	100 mA @ 28 VDC
Warranty:	2-year warranty against defects in materials and workmanship from date of purchase



All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.



Specifications

Height:

7.5"

Width:

2.6" (tapering to 1.9")

Thickness:

1.7"

Weight:

18.2 oz

Cord Length:

8'

Color:

White

Case Material:

UL 94V-0 Flame Rated, Latex Free ABS

Speaker Material:

Moisture-resistant mylar cone

Electrical Rating:

24 VDC at 100 mA

Standards:

Tested to UL 1069 O2 SafeTesting with oxygen enriched environment marking
RoHS Compliant
FDA Listed

Warranty:

2-year warranty against defects in materials and workmanship from date of purchase

LATEX-FREE

OXYGEN SAFE

QUANTITY DISCOUNTS
CALL FOR PRICING


2
YEAR WARRANTY

Slimline™

Crest Slimline™ Pillow Speakers

Slimline Pillow Speakers are designed for each hospital's specific needs and fit many nurse call systems. These units feature a slim “easy to hold” design, with many options available. Slimline Pillow Speakers are priced far below original equipment cost and are designed to make repairs quick and easy.

- Slim, easy to hold design
- Comes standard with 8' cord length



Part No.	System	Nurse Call	TV Control	Other Features	Price
PSS38010801†	Jeron Pro-Alert 640 and Provider 680		Yes		82.57
PSS38020802†	Jeron Pro-Alert 640 and Provider 680	Nurse	Yes		82.57
PSS38040801†	Jeron Pro-Alert 640 and Provider 680	Nurse	Yes	1 light	100.54

Replacement Cable Assemblies for Crest Slimline™

Replacement cable assemblies for Slimline are easy to replace. Standard white cables are 8' with molded-in strain relief. To verify the correct cable for your pillow speaker or for replacement assistance, please call Crest Technical Support.


Slimline Cable Assemblies					
Part No.	Style	System	Color	Length	Price
PS3801W-08	8-pin DIN	Jeron Pro-Alert 640 and Provider 680	White	8'	32.13




Crest EconoCall™ Call Cords

The EconoCall delivers a clean, basic design with reliable function at a price any facility can afford.


- Our most economically priced call cord
- Repairable button for fast fixes
- Plug styles available to fit most systems
- RoHS 2 Compliant



Part No.	System	Style	Color	Length	From
117700-__†	Any station with ¼" jack or Type B or Y adapter	Single	White	7' or 12'	18.09



Part No.	System	Style	Color	Length	From
7700W-__†	Any station with ¼" jack or Type B or Y adapter	Single	White	7', 10', 12', 15' or 20'	14.58
7700B-__†	Any station with ¼" jack or Type B or Y adapter	Single	Beige	7', 10', or 12'	18.92
7700G-__†	Any station with ¼" jack or Type B or Y adapter	Single	Grey	7', 10' or 12'	18.83
7700DW-__†	Any station with ¼" jack or Type B or Y adapter	Double	White	7', 10', 12', 15' or 20'	26.37
7700DB-__†	Any station with ¼" jack or Type B or Y adapter	Double	Beige	10' or 12'	31.54
7700DG-__†	Any station with ¼" jack or Type B or Y adapter	Double	Grey	10' or 12'	31.54



Part No.	System	Style	Color	Length	Price
117738†	Jeron Pro-Alert 640 and Provider 680	Single	White	7'	23.95
117738-12†	Jeron Pro-Alert 640 and Provider 680	Single	White	12'	38.53

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.



Crest CleanCall™ Call Cords

Maintaining a clean environment can be difficult, but CleanCall simplifies the job. The completely sealed membrane style button can be easily wiped down to prevent the spread of germs, making it a great option for infection control in any facility.



- Smooth membrane is easy to sanitize
- Easy grasp winged design fits comfortably in the hand
- Highly visible, high-contrast button with nurse call symbol

LATEX-FREE

OXYGEN SAFE



1 YEAR WARRANTY

¼" Phone Plug
2-conductor



Part No.	System	Style	Color	Length	From
8900W-__ [†]	Any station with ¼" jack or Type B or Y adapter	Single	White	8', 10', 12' or 15'	27.26
8900DW-10 [†]	Any station with ¼" jack or Type B or Y adapter	Double	White	10'	61.97

8-pin DIN Plug



Part No.	System	Style	Color	Length	Price
8938W-8 [†]	Jeron Pro-Alert 640 and Provider 680	Single	White	8'	34.54

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.



Crest DuraCall™ Call Cords

When durability and feature versatility are a must, look no further than the DuraCall Call Cord line. With its rugged design, DuraCall has survived drops, submersion in water and even being driven over by a car.

- Easy-grip winged pendant
- Tactile "click" raised button
- Easy to clean
- Wide variety of plug styles & call cord features available
- Short 6" length available for safety in anti-ligature environments



RECOGNIZED COMPONENT

LATEX-FREE

OXYGEN SAFE

2 YEAR WARRANTY



¼" Phone Plug
2-conductor





Part No.	System	Style	Color	Length	From
9900W-5 [†]	Any station with ¼" jack or Type B or Y adapter	Anti-Ligature	White	6"	33.82
9900W-__ [†]	Any station with ¼" jack or Type B or Y adapter	Single	White	7', 10', 12', 15' or 20'	16.05
9900DW-__ [†]	Any station with ¼" jack or Type B or Y adapter	Double	White	7', 10', 12', 15' or 20'	26.95
9900CLW-16 [†]	Any station with ¼" jack or Type B or Y adapter	Single, Coiled	White	16'	29.80

8-pin DIN Plug



Part No.	System	Style	Color	Length	Price
9938W-7 [†]	Jeron Pro-Alert 640 and Provider 680	Single	White	7'	27.62
9938W-10 [†]	Jeron Pro-Alert 640 and Provider 680	Single	White	10'	30.96
9938W-15 [†]	Jeron Pro-Alert 640 and Provider 680	Single	White	15'	42.06



All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.

Crest DuraColor™ Call Cords

Crest’s neon green pendant in the DuraCall Call Cord line is the latest option for nurse call visibility and ease of use. Bright and vibrant, the pendant is ideal for those with poor eyesight or for those who may have difficulty distinguishing it from bedding.

- High contrast neon color for easy identification
- Prominent color assists patients & staff
- White cable

1/4" Phone Plug
2-conductor







Part No.	System	Other Features	Style	Color	Length	From
9900N-__ †	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	High-Contrast Color Pendant	Single	White	7' or 10'	22.39

Crest DuraGlo™ Call Cords

DuraGlo is the only call cord in the industry with glowing bands to help your residents locate their call cord in dark rooms! DuraGlo can glow up to four hours in total darkness. Built on the same solid platform as our popular DuraCall call cord. Works with any nurse call system.

- Easily locate in dark rooms
- Recharges with light exposure
- Requires no power or batteries

1/4" Phone Plug
2-conductor







Part No.	System	Other Features	Style	Color	Length	From
9900L-__ †	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Glow-in-the-Dark Pendant	Single	White	7' or 10'	33.29

Crest DuraLight™ Call Cords

DuraLight Call cords for 1/4" phone plug systems have LEDs that momentarily light up when the button is depressed (may not work with every system).

- Bright red LEDs light up the pendant
- Requires no batteries
- FDA Listed

1/4" Phone Plug
2-conductor





Part No.	System	Other Features	Style	Color	Length	From
9900R-__ †	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Red LEDs	Single	White	7' or 10'	32.52

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.



Crest SoftCall™ Momentary Call Cords



SoftCall provides a durable and affordable option for nurse call and light controls. The ergonomic pendant fits comfortably in the hand, and its silicone skin provides a reassuring, “soft” feel. The waterproof pendant can easily be wiped down for infection control. The high-contrast color makes the SoftCall pendant easy to find in bedding.

- Soft, silicone pendant
- Easy to clean
- Durable strain relief
- IP47 sealed & waterproof

LATEX-FREE



2
YEAR
WARRANTY

1/4" Phone Plug
2-conductor



Part No.	System	Style	Color	Length	Price
SC44B-01 †	Any station with 1/4" jack or Type B or Y adapter	Nurse	Blue	8'	43.55
SC44B-02 †	Any station with 1/4" jack or Type B or Y adapter	Nurse	Blue	12'	45.42

8-pin DIN Plug




Part No.	System	Style	Color	Length	Price
SC38B-01 †	Jeron Pro-Alert 640 and Provider 680	Nurse	Blue	8'	60.80
SC38C-01 †	Jeron Pro-Alert 640 and Provider 680	Nurse, 1 Light	Blue	8'	62.09

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.

BreathCall® Air Activated Call Cords

BreathCall is for patients with little to no motor skills and is activated by simply breathing into the disposable straw and filter assembly. It has a flexible gooseneck with universal mounting clamp. The Mini BreathCall offers the same features of the standard with a more streamlined design to fit smaller areas. All BreathCalls come with a 9' cord, clamp and filter and straw kit. Standard BreathCalls have 36" black tubing from clamp to mouthpiece, and Mini BreathCalls have 28" white tubing. For single patient use.

- Flexible conduit bends to almost any position
- Sturdy clamp easily attaches to headboards
- Locking or non-locking switch closure available
- Dual controls for nurse, TV or radio available (11201-09D)



Part No.	System	Style	Color	Length	Price
11201-09	Any station with ¼" jack or Type B or Y adapter	Single, Non-Locking	Black	9'	167.30
11201-09D	Any station with ¼" jack or Type B or Y adapter	Double, Non-Locking	Black	9'	220.86
16001-09	Any station with ¼" jack or Type B or Y adapter	Mini, Single, Non-Locking	White	9'	133.67




11201-09


Crest GeriCall™ Call Cords

The GeriCall Geriatric Call Cord is a round, flattened momentary call cord pad that works great for residents who have limited hand dexterity. This call cord utilizes a mechanical switch in the pad that can be activated by slight pressure from the hand, arm or body, or even can be positioned under the resident's chin.

- Large pad for easy activation
- Includes security clip
- 3" diameter pad



Part No.	System	Style	Color	Length	From
118844C[†]	Any station with ¼" jack or Type B or Y adapter	Single	Grey	8'	34.48
118844C-__[†]	Any station with ¼" jack or Type B or Y adapter	Single	Grey	10' or 15'	37.45
118844D-10[†]	Any station with ¼" jack or Type B or Y adapter	Double	Grey	10'	88.26



Part No.	System	Style	Color	Length	Price
118838C-10[†]	Jeron Pro-Alert 640 and Provider 680	Single	Grey	10'	51.97




All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.


PressCall® Pneumatic Call Cords

PressCall is a pneumatic air-activated call cord designed for patients who have difficulty activating standard call cords. The PressCall has an oval shaped, soft vinyl bulb design that fits easily in the hand.

- For patients with sufficient motor skills too press bulb
- Great for patients unable to operate regular call cords
- Easy to squeeze bulb style



Part No.	System	Style	Color	Length	From
11372__	Any station with ¼" jack or Type B or Y adapter	Single, Non-Locking	White or Grey	6', 8', 10', 12' or 15'	42.65
1159__	Any station with ¼" jack or Type B or Y adapter	Double, Non-Locking	White or Grey	6', 8', or 12'	60.85
11373__	Any station with ¼" jack or Type B or Y adapter	Double, Non-Locking	White or Grey	8' or 10'	67.58




Part No.	System	Style	Color	Length	Price
113776	Jeron Pro-Alert 640 and Provider 680	Single, Non-Locking	Grey	8'	77.60


PadCall® Pneumatic Call Cords

PadCall is a pneumatic air-activated call cord designed for patients who have difficulty activating standard call cords. They have a large circular, soft vinyl pad with gripping feet that prevent it from slipping off any surface.

- For patients with sufficient motor skills too press or lean on pad
- Great for patients unable to operate regular call cords
- Stay-put 3 ½" wide pad



Part No.	System	Style	Color	Length	From
1137__	Any station with ¼" jack or Type B or Y adapter	Single, Non-Locking	White or Grey	6', 8', 10', 12' or 15'	From
1159__	Any station with ¼" jack or Type B or Y adapter	Double, Non-Locking	White or Grey	6', 8', or 12'	83.45
1137__	Any station with ¼" jack or Type B or Y adapter	Double, Non-Locking	White or Grey	8' or 10'	102.25



Part No.	System	Style	Color	Length	Price
113777	Jeron Pro-Alert 640 and Provider 680	Single, Non-Locking	Grey	8'	88.71

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors V1. Revised 07/25/23.



Crest TouchCall™ Call Cords

TouchCall provides reliable and versatile nurse call actuation for patients with limited dexterity. The TouchCall features a large, ultra-sensitive touch plate surface that can be easily pressed by hand, foot, elbow, chin, etc. The tent design with clip and Velcro® attachment methods allow for multiple positioning possibilities to make it simple and convenient for any patient to use.

- Standard 10' cord length to suit a variety of applications
- Classic grey color for easy visibility

LATEX-FREE

OXYGEN SAFE

2 YEAR WARRANTY

1/4" Phone Plug
2-conductor

Part No.	System	Other Features	Style	Color	Length	Price
113673-4410 [†]	Any station with 1/4" jack or Type B or Y adapter	Ultra-Sensitive Touch Plate	Single	Grey	10'	64.44

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.



Adapters

- Type B Adapters** convert your existing receptacle to the same receptacle and 1/4" phone receptacle. Include 9" extension cord and plug.
- Type Y Adapters** convert your existing receptacle to two 1/4" phone receptacles.
- Type Y2 Adapters** convert your existing receptacle to one female receptacle of the same style and one 1/4" phone receptacle.
- Type Y3 Adapters** convert your existing receptacle to two of your existing receptacles.
- Type D Adapters** convert your existing receptacle to one 1/4" phone receptacle.

Part No.	System	Style	Description	Price
9520	Jeron Pro-Alert 640 and Provider 680	Type B	8-pin DIN to 8-pin DIN and 1/4" 2C	93.68
CA38X00N2-01	Jeron Pro-Alert 640 and Provider 680	Type Y	8-pin DIN to (2) 1/4" 2C	66.87
CA38S38N1-01	Jeron Pro-Alert 640 and Provider 680	Type Y2	8-pin DIN to 8-pin DIN and 1/4" 2C	59.40
CA38D38N1-01	Jeron Pro-Alert 640 and Provider 680	Type Y3	8-pin DIN to (2) 8-pin DIN	33.94
114534	Jeron Pro-Alert 640 and Provider 680	Type D	8-pin DIN to 1/4" 2C	14.32



Station Savers

Station Savers create a break-away point between cord and station to help prevent damage. Disconnection of cord triggers cord-out on stations with this feature.

Part No.	System	Description	Length	Price
108090-01	Any station with 1/4" jack or Type B or Y adapter	1/4", Right Angle	1'	38.50
107397	Jeron Pro-Alert 640 and Provider 680	8-pin DIN	1'	25.03
114711-01	Jeron Pro-Alert 640 and Provider 680	8-pin DIN, Right Angle	1'	25.27

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V1. Revised 07/25/23.

Call Cord Extenders

Call Cord Extenders allow you to lengthen any 2-conductor, ¼" phone plug nurse call cord.

- Disconnection of cord triggers cord-out on stations with this feature
- Multiple lengths to fit any application



114711-10

Part No.	System	Color	Length	Price
108090	Any station with ¼" jack or Type B or Y adapter	White	8'	24.03
108090-12	Any station with ¼" jack or Type B or Y adapter	White	12'	24.55
114711-10	Jeron Pro-Alert 640 and Provider 680	White	10'	28.48
107397-10	Jeron Pro-Alert 640 and Provider 680	White	10'	39.95

Dummy Plugs

Dummy plugs provide proper termination too eliminate cord-out calls from being placed when a call cord, pillow speaker or bed is disconnected from the wall station.



Part No.	System	Price
107248D	Jeron Pro-Alert 640 and Provider 680	24.03



B-32877



B-32877P



B-33509

Replacement Receptacles

Part No.	System	Description	Price
B-32877	Jeron Pro-Alert 640 and Provider 680	8-pin DIN Receptacle, panel mount	8.87
B-32877P	Jeron Pro-Alert 640 and Provider 680	8-pin DIN Receptacle, on 1-gang stainless steel plate	49.65
B-33509	Jeron Pro-Alert 640 and Provider 680	8-pin DIN, Receptacle, PC mount	11.87

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V1. Revised 07/25/23.



6420C

6422C



B-31993 Red



E-11657 Green LED

Patient Stations

Crest replacement stations for Jeron® Pro-Alert™ 640 Nurse Call System provide two-way communication between staff and patients.

Part No.	System	Description	Price
6420C [†]	Jeron Pro-Alert 640 and Provider 680	Patient Station, single jack	208.40
6422C [†]	Jeron Pro-Alert 640 and Provider 680	Patient Station, dual jack	255.39



B-31993

Waterproof Pullcord Bath Station

Crest replacement Pullcord Bath Station for Jeron® Pro-Alert™ 640 system. For use in bathrooms or other wet environments. Mounts on a 1-gang box, with a plastic faceplate.



Part No.	System	Description	Price
6432C [†]	Jeron Pro-Alert 640 and Provider 680	Pullcord Bath Station, waterproof	61.75

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 07/25/23.