



Crest Healthcare

Jeron® Provider® 790 Part List

Your guide to Jeron® Provider 790 system parts available from Crest that are durable, economical and designed with infection control in mind.

CREST
HEALTHCARE®

© 2023 Crest Healthcare



Enhancing the environment of care for those in need and those who care

Crest Healthcare is a leading provider of a wide range of healthcare products and technical services. Historically specialized in high-quality nurse call replacement parts, and today producing nurse call system components for nearly every current system, Crest has been a trusted source for innovative product development and production since 1967.

With customer satisfaction as our driving goal, Crest products provide quality, value, innovation, variety and durability. Working together to guarantee a positive experience, our in-house team includes customer service, technical support, engineering, manufacturing and quality assurance experts.

Engineering

Our engineering team is continually creating and testing new communication products. With the use of destructive and non-destructive in-house testing and prototype evaluations, Crest ensures every product provides durability and a long service life. Custom services are available to OEMs and other commercial customers.

Manufacturing

We build and assemble all the Crest products in-house and work closely with our manufacturing partners to ensure our same high quality standards are met.


Quality Assurance

We conduct testing on 100% of finished products we manufacture prior to shipment. Crest products meet or exceed demanding safety standards up to and including UL 1069 Recognition.

Quality Certifications

Crest strives to offer you high quality products to support your facility needs. Crest Healthcare is ISO 13485:2016 certified provider of Pillow Speakers, Nurse Call Stations, Nurse Call Cords, and Bed Control Equipment. Our certification represents our commitment to on-going quality improvements and a dedication to higher levels of performance, service and quality products. We received our first certification in 1999 (ISO 9001) and Crest continues to improve quality in all aspects of designing, manufacturing and customer service.



 800.328.8908

 800.369.9207

 CrestHealthcare.com

Jeron® Provider® 790

Nurse Call Components

Pillow Speakers

Elite Digital™	4
Elite Infrared™	5
Replacement Cable Assemblies	5

Call Cords

Momentary Call Cords

DuraCall™	6
PermaLight™	6

Specialty Call Cords

GeriCall™	7
-----------------	---




EliteDigital™

Crest Backlit Elite Digital™ Pillow Speakers

Designed for use with digital TVs, Autosense™ Technology automatically identifies the brand of digital healthcare TV and sends the correct signal. They can be moved to any room, regardless of healthcare television brand.

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille



Specifications

Height:

9"

Width:

4 3⁄5"

Thickness:

1 3⁄4"

Cord Length:

8'

Case Material:

ABS, UL94V-0

Cord Connector:

Plug-and-Play (solderless)

Power Requirements:

Two AAA alkaline batteries required (included) (Infrared only)

Switch Rating:

30 VDC @ 100 mA

Warranty:

2-year warranty against defects in materials and workmanship from date of purchase

LATEX-FREE

OXYGEN SAFE

PATENTED

QUANTITY DISCOUNTS
CALL FOR PRICING

2
YEAR WARRANTY

ETL
US

Part No.	System	Nurse Call	TV Control	Other Features	Price
XN5301XPBBXA14N1†	Jeron Provider 790	Nurse, pain, toilet, water	10-key	2 lights, headphone, backlit	171.99

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors.† Item is FDA listed. V1. Revised 05/10/23.

EliteInfrared™

Crest Backlit Elite Infrared™ Pillow Speakers

The Crest Elite Infrared Pillow Speaker is the only pillow speaker that combines nurse call and TV with external DVD or Set Top Box (STB) controls. It combines the ease of a universal remote with traditional pillow speaker functionality such as nurse call and environmental controls. It is also useful where there is no cross-room wiring or with TVs that do not have a pillow speaker port. Uses two AAA batteries (included).

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille

8-pin DIN Plug



Specifications

Height:	9"
Width:	4 3/5"
Thickness:	1 3/4"
Cord Length:	8'
Case Material:	ABS, UL94V-0
Cord Connector:	Plug-and-Play (solderless)
Power Requirements:	Two AAA alkaline batteries required (included) (Infrared only)
Switch Rating:	30 VDC @ 100 mA
Warranty:	2-year warranty against defects in materials and workmanship from date of purchase

LATEX-FREE

OXYGEN SAFE



PATENTED

QUANTITY DISCOUNTS
CALL FOR PRICING

2
YEAR
WARRANTY



Part No.	System	Nurse Call	TV Control	Other Features	Price
XT5301PPBBXXZ14N1†	Jeron Provider 790	Nurse, pain, toilet, water	10-key	2 lights, headphone, backlit	193.16

CONFIGURATION QUESTIONS?

CONTACT TECH SUPPORT

1-833-421-0294

Replacement Cable Assemblies for Elite Pillow Speakers

Replacement pillow speaker cable assemblies for Elite Digital™ and Elite Infrared™ are easy to replace using "plug-and-play" connectors. Standard white cables are 8' with molded-in strain relief. To verify the correct cable for your pillow speaker or for replacement assistance, please contact Crest Technical Support.

8-pin DIN Plug



Part No.	System	Color	Length	Price
XN5301W-08	Jeron Provider 790	White	8'	19.58

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 05/10/23.

Crest DuraCall™ & PermaLight™ Momentary Call Cords




Crest DuraCall™ Call Cords

When durability and feature versatility are a must, look no further than the DuraCall Call Cord line. With its rugged design, DuraCall has survived drops, submersion in water and even being driven over by a car.

- Easy-grip winged pendant
- Tactile "click" raised button
- Easy to clean
- Wide variety of plug styles & call cord features available
- Short 6" length available for safety in anti-ligature environments



Durable Waterproof Design





Part No.	System	Style	Color	Length	Price
9953W-7†	Jeron Provider 790	Single	White	7'	27.59
9953W-10†	Jeron Provider 790	Single	White	10'	28.51

Crest PermaLight™ Call Cords

The Jeron Provider 790 PermaLight has LEDs that are continually illuminated for visibility in low-light environments. Built on the same solid platform as our popular DuraCall Call Cords.

- Bright red LEDs light up the pendant
- Requires no batteries
- Easily locate in dark rooms



Part No.	System	Other Features	Style	Color	Length	Price
9953P-7†	Jeron Provider 790	Red LEDs	Single	White	7'	33.58

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 05/10/23.

Crest GeriCall™ Call Cords

The GeriCall Geriatric Call Cord is a round, flattened momentary call cord pad that works great for residents who have limited hand dexterity. This call cord utilizes a mechanical switch in the pad that can be activated by slight pressure from the hand, arm or body, or even can be positioned under the resident's chin.

- Large pad for easy activation
- Includes security clip
- 3" diameter pad



8-pin DIN Plug



LATEX-FREE

OXYGEN
SAFE

2
YEAR
WARRANTY

Part No.	System	Style	Color	Length	Price
118853-10†	Jeron Provider 790	Single	Grey	10'	54.12



www.CrestHealthcare.com

FULL CATALOG AVAILABLE ONLINE



- Over 3,000 products available
- Order all your facility products
- Bookmark, print or email pages

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V1. Revised 05/10/23.