



Crest Healthcare

Rauland® Responder® III+, IV, & 4000 Parts List

Your guide to Rauland® Responder® III+, IV and 4000 system parts available from Crest that are durable, economical and designed with infection control in mind.

CREST
HEALTHCARE®

© 2023 Crest Healthcare



Enhancing the environment of care for those in need and those who care

Crest Healthcare is a leading provider of a wide range of healthcare products and technical services. Historically specialized in high-quality nurse call replacement parts, and today producing nurse call system components for nearly every current system, Crest has been a trusted source for innovative product development and production since 1967.

With customer satisfaction as our driving goal, Crest products provide quality, value, innovation, variety and durability. Working together to guarantee a positive experience, our in-house team includes customer service, technical support, engineering, manufacturing and quality assurance experts.

Engineering

Our engineering team is continually creating and testing new communication products. With the use of destructive and non-destructive in-house testing and prototype evaluations, Crest ensures every product provides durability and a long service life. Custom services are available to OEMs and other commercial customers.

Manufacturing

We build and assemble all the Crest products in-house and work closely with our manufacturing partners to ensure our same high quality standards are met.

Quality Assurance

We conduct testing on 100% of finished products we manufacture prior to shipment. Crest products meet or exceed demanding safety standards up to and including UL 1069 Recognition.

Quality Certifications

Crest strives to offer you high quality products to support your facility needs. Crest Healthcare is ISO 13485:2016 certified provider of Pillow Speakers, Nurse Call Stations, Nurse Call Cords, and Bed Control Equipment. Our certification represents our commitment to on-going quality improvements and a dedication to higher levels of performance, service and quality products. We received our first certification in 1999 (ISO 9001) and Crest continues to improve quality in all aspects of designing, manufacturing and customer service.



 800.328.8908

 800.369.9207

 CrestHealthcare.com

Rauland® Responder® III+, IV & 4000
Nurse Call Components

Pillow Speakers

Elite Digital™ 4 - 5

Elite Infrared™ 4 - 5

Replacement Cable Assemblies 6

Pillow Speaker Tester 6

Remote Call-IR™ 7

Slimline™ 8

Call Cords

Momentary Call Cords

EconoCall™ 9

CleanCall™ 10

DuraCall™ 11

DuraColor™ 12

DuraGlo™ 12

DuraLight™ 12

SoftCall™ 13

JumboCall™ 13

Specialty Call Cords

BreathCall® 14

GeriCall™ 15

PressCall® 16

PadCall® 16

TouchCall™ 17

Additional Components

Magnetic Break-A-Way Cable System 18

Call Cord Extensions 19

Station Savers..... 19

Adapters 19

Replacement Receptacles 20

Corridor Lights..... 20

Rauland Responder IV Stations 21

Rauland Responder 4000 Stations..... 22

Rauland Responder III+ Stations 23



EliteDigital™

Crest Elite Digital™ Pillow Speakers

Designed for use with healthcare TVs. Autosense™ Technology automatically identifies the brand of healthcare TV and sends the correct signal, eliminating programming hassle. They can be moved to any room, regardless of healthcare television brand.

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille

EliteInfrared™

Crest Elite Infrared™ Pillow Speakers

The Crest Elite Infrared Pillow Speaker is the only pillow speaker that combines nurse call and TV with external DVD or Set Top Box (STB) controls. It combines the ease of a universal remote with traditional pillow speaker functionality such as nurse call and environmental controls. It is also useful where there is no cross-room wiring or with TVs that do not have a pillow speaker port. Uses two AAA batteries (included).

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille

Specifications

Height:

9"

Width:

4 3⁄8"

Thickness:

1 3⁄4"

Cord Length:

8'

Case Material:

ABS, UL94V-0

Cord Connector:

Plug-and-Play (solderless)

Power Requirements:

Two AAA alkaline batteries required (included) (Infrared only)

Switch Rating:

30 VDC @ 100 mA

Warranty:

2-year warranty against defects in materials and workmanship from date of purchase

LATEX-FREE

OXYGEN SAFE

PATENTED

QUANTITY DISCOUNTS
CALL FOR PRICING

2
YEAR WARRANTY

ETL
US



All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V20. Revised 04/24/23.

EliteDigital™

8-pin DIN Plug

Part No.	System	Nurse Call	TV Control	Other Features	Price
XN0401XBA014A10N1†	Rauland Responder III+, IV, and 4000	Nurse	10-key		119.95
XN0401XCA014A10N1†	Rauland Responder III+, IV, and 4000	Nurse	10-key	1 light	121.35
XN0401XEA014A10N1†	Rauland Responder III+, IV, and 4000	Nurse	10-key	2 lights	119.95
XN0401EBA014A10N1†	Rauland Responder III+, IV, and 4000	Nurse	Up/Down		115.45
XN0401ECA014A10N1†	Rauland Responder III+, IV, and 4000	Nurse	Up/Down	1 light	115.45
XN0401EEA014A10N1†	Rauland Responder III+, IV, and 4000	Nurse	Up/Down	2 lights	115.45

1⁄4" Phone Plug
2-conductor

EliteInfrared™

Part No.	System	Nurse Call	TV Control	Other Features	Price
XT4401XBA600Z10N1†	Any station with 1⁄4" jack or Type D, B, Y, or Y2 adapter	Nurse	10-key		148.46
XT4401DBA600Z10N1†	Any station with 1⁄4" jack or Type D, B, Y, or Y2 adapter	Nurse	10-key	DVD	172.03

8-pin DIN Plug

Part No.	System	Nurse Call	TV Control	Other Features	Price
XT0401XBA014Z10N1†	Rauland Responder III+, IV, and 4000	Nurse	10-key		148.46
XT0401XCA014Z10N1†	Rauland Responder III+, IV, and 4000	Nurse	10-key	1 light	148.46
XT0401XEA014Z10N1†	Rauland Responder III+, IV, and 4000	Nurse	10-key	2 lights	148.46
XT0401EBA014Z10N1†	Rauland Responder III+, IV, and 4000	Nurse	Up/Down		142.57
XT0401ECA014Z10N1†	Rauland Responder III+, IV, and 4000	Nurse	Up/Down	1 light	142.57
XT0401EEA014Z10N1†	Rauland Responder III+, IV, and 4000	Nurse	Up/Down	2 lights	142.57

Questions? Contact Us For Assistance!

Contact Crest Technical Support for further assistance.

1-833-421-0294
techsupport@cresthealthcare.com



All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.





XN4401W-08

Replacement Cable Assemblies for Elite Pillow Speakers



Replacement pillow speaker cable assemblies for Elite Digital™ and Elite Infrared™ are easy to replace using simple "plug-and-play" connectors. Standard white cables are 8' with molded-in strain relief. To verify the correct cable for your pillow speaker or for replacement assistance, please call Crest Technical Support at 1-833-421-0294.

1/4" Phone Plug
2-conductor




Part No.	System	Color	Length	Price
XN4401W-08	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	White	8'	25.95

8-pin DIN Plug



Part No.	System	Color	Length	Price
XN0401W-08	Rauland Responder III+, IV, and 4000	White	8'	32.13



Bottom View

Crest Pillow Speaker Tester

Test functionality of analog, digital and infrared pillow speakers as well as bed controls and call cords with this compact and portable tester. Please call Crest's Technical Support team with your pillow speaker model number for assistance in obtaining a module to fit your system.

- Operates with 12 VDC or 9 V battery (both included)
- Tester works with bed controls & call cords
- 37 Pin Module connection to pillow speaker tester
- Modules to fit nurse call systems sold separately

Specifications

Length:	5 3/4"
Width:	4 3/5"
Height:	1 3/4"
Power Requirements:	9 V Battery or 12 VDC adapter (both included)



Part No.	System	Description	Price
115642	Rauland Responder III and 3000	Speaker Tester, Not Compatible with Rauland Responder 5	268.21

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V20. Revised 04/24/23.





Crest Remote Call-IR™ TV Control & Nurse Call

Replace your call cord and TV remote control with one simple device. Remote Call-IR controls the TV using IR signals like a universal remote while still providing the safety of reliable wired nurse call. Remote controls are easily lost or broken, but the Remote Call-IR plugs into your nurse call station to ensure it stays within the room and is durable enough to withstand drops and daily use. The 10-key overlay style is great for emergency rooms or private patient rooms, where the room is not wired for TV sound or the sound control of a pillow speaker is not needed. The large button overlay is perfect for kids, geriatric patients or residents with limited dexterity. Features simplified control of TV power, channel and volume with big, easy-to-see buttons. Both styles feature an on-board library of codes that can be used to program it to work with almost any TV. Uses two AA batteries (not included) that last as long as in typical TV remotes.

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille
- Bed A & Bed B functionality
- Fits Crest Pillow Speaker Keeper

1/4" Phone Plug
2-conductor




Part No.	System	Nurse Call	Style	Length	Price
XP441E0B600XXX2 †	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Nurse	Large Button	10'	81.95
XP441X0B600XXX2 †	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Nurse	10-key	10'	89.95

Replacement Cable for Remote Call IR

Part No.	System	Length	Price
XP4401W-10	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	10'	33.91

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.

 Programming instructions and code list can be found at www.cresthealthcare.com under the "Resources" tab

Specifications

Height:	6 1/2"
Width:	Tapers, 2 3/8" - 2 1/4"
Thickness:	1 3/8" at thickest point
Weight:	13 7/10 oz (with batteries)
Cord Length:	10'
Case Material:	ABS, Flame rated (UL94V-0)
Cord & Strain Relief:	Flexible molded PVC
Cord Connector:	Plug-and-Play (solderless)
Plug Style:	2-conductor or 3-conductor 1/4" phone plug (call Crest for available adapters)
Power Requirements:	Two AA alkaline batteries (sold separately)
Switch Rating:	100 mA @ 28 VDC
Warranty:	2-year warranty against defects in materials and workmanship from date of purchase

LATEX-FREE

OXYGEN SAFE

QUANTITY DISCOUNTS
CALL FOR PRICING

2 YEAR WARRANTY

114576




Slimline™

Crest Slimline™ Pillow Speakers

The Slimline Pillow Speaker is designed for each hospital's specific needs and can fit many nurse call systems. These units feature a slim "easy to hold" design, with many options available. Slimline pillow speakers are priced far below original equipment cost and are designed to make repairs quick and easy.


- Slim, easy to hold design
- Comes standard with 8' cord length



Part No.	System	Nurse Call	TV Control	Other Features	Price
PSS04010801 †	Rauland Responder III+, IV, and 4000		Yes		82.57
PSS04020801 †	Rauland Responder III+, IV, and 4000	Nurse	Yes		82.57
PSS04040801 †	Rauland Responder III+, IV, and 4000	Nurse	Yes	1 light	89.05

Crest Slimline™ Cable Assemblies

Replacement cable assemblies for Slimline are easy to replace. Standard white cables are 8' with molded-in strain relief. To verify the correct cable for your pillow speaker or for replacement assistance, please call Crest Technical Support at 1-833-421-0294.



Part No.	System	Color	Length	Price
PS0401W-08	Rauland Responder III+, IV, and 4000	White	8'	32.13

Specifications

Height:	7.5"
Width:	2.6" (tapering to 1.9")
Thickness:	1.7"
Weight:	18.2 oz
Cord Length:	8'
Color:	White
Case Material:	UL 94V-0 Flame Rated, Latex Free ABS
Speaker Material:	Moisture-resistant mylar cone
Electrical Rating:	24 VDC at 100 mA
Standards:	Tested to UL 1069 O2 SafeTesting with oxygen enriched environment marking RoHS Compliant FDA Listed
Warranty:	2-year warranty against defects in materials and workmanship from date of purchase



Crest EconoCall™ Call Cords

If you're looking for an affordable and reliable call cord, the EconoCall delivers. The EconoCall combines a clean, basic design with reliable function at a price any facility can afford. The smooth pendant shape fits comfortably in the hand.

- Our most economically priced call cord
- Repairable button for fast fixes
- Plug styles available to fit most systems
- RoHS 2 Compliant




Straight
¼" Phone Plug
2-conductor

Part No.	System	Style	Color	Length	From
117700-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	White	7' or 12'	18.09

¼" Phone Plug
2-conductor

Part No.	System	Style	Color	Length	From
7700W-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	White	7', 10', 12', 15' or 20'	14.58
7700B-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	Beige	7', 10', or 12'	15.80
7700G-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	Grey	7', 10' or 12'	18.83
7700DW-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Double	White	7', 10', 12', 15' or 20'	26.37
7700DB-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Double	Beige	10' or 12'	31.54
7700DG-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Double	Grey	10' or 12'	31.54



Part No.	System	Style	Color	Length	From
117704-__ †	Rauland Responder III+, IV, and 4000	Single	White	7' or 12'	23.95
117704D-7 †	Rauland Responder III+, IV, and 4000	Double	White	7'	47.88

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.



Crest CleanCall™ Call Cords

Maintaining a clean environment can be difficult, but the CleanCall helps simplify the job. The completely sealed, membrane style button can be easily wiped down to help prevent the spread of germs, making it a great option for infection control in any facility.

- Smooth membrane is easier to sanitize than regular buttons
- Easy-grasp winged design fits comfortably in the hand
- High-contrast red button with nurse call symbol is highly visible
- Tactile "click" button



LATEX-FREE

OXYGEN SAFE

1 YEAR WARRANTY





¼" Phone Plug
2-conductor



Part No.	System	Style	Color	Length	From
8900W-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	White	8', 10', 12' or 15'	26.09
8900DW-10 †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Double	White	10'	59.39

8-pin DIN Plug



Part No.	System	Style	Color	Length	Price
8904W-8 †	Rauland Responder III+, IV, and 4000	Single	White	8'	32.86

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.



Crest DuraCall™ Call Cords



When durability and feature versatility are a must, look no further than the DuraCall Call Cord line. With its rugged design, DuraCall has survived drops, submersion in water and even being driven over by a car.

- Easy-grip winged pendant
- Raised button
- Easy to clean & sanitize
- Wide variety of plug styles & call cord features available
- Short 6" length available for safety in anti ligature environments





¼" Phone Plug
2-conductor



Part No.	System	Style	Color	Length	From
9900W-5 †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Anti-Ligature	White	6"	26.00
9900W-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	White	7', 10', 12', 15' or 20'	16.05
9900DW-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Double	White	7', 10', 12', 15' or 20'	26.95
9900CLW-16 †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single, Coiled	White	16'	29.80

8-pin DIN Plug



Part No.	System	Style	Color	Length	From
9904W-__ †	Rauland Responder III+, IV, and 4000	Single	White	7' or 10'	25.45
9904DW-7 †	Rauland Responder III+, IV and 4000	Double	White	7'	40.71


All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.

Crest DuraColor™, DuraGlo™ & DuraLight™ Momentary Call Cords

Crest DuraColor™ Call Cords

DuraColor features a bright, neon green pendant that is ideal for those with poor eye-sight or for those who may have difficulty distinguishing it from bedding.

- High contrast neon color for easy identification
- Prominent color assists patients & staff
- White cable




Part No.	System	Other Features	Style	Color	Length	From
9900N-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	High-Contrast Color Pendant	Single	White	7' or 10'	22.39

Crest DuraGlo™ Call Cords

DuraGlo is the only call cord in the industry with glowing bands to help your residents locate their call cord in dark rooms! DuraGlo can glow up to four hours in total darkness. Built on the same solid platform as our popular DuraCall Call Cords. Works with any nurse call system.

- Easily locate in dark rooms
- Recharges with light exposure
- Requires no power or batteries




Part No.	System	Other Features	Style	Color	Length	From
9900L-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Glow-in-the-Dark Pendant	Single	White	7' or 10'	33.29

Crest DuraLight™ Call Cords

DuraLight Call cords for ¼" phone plug systems have LEDs that momentarily light up when the button is depressed (may not work with every system).

- Bright red LEDs light up the pendant
- Requires no batteries
- FDA Listed



Part No.	System	Other Features	Style	Color	Length	From
9900R-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Red LEDs	Single	White	7' or 10'	32.52


All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.

Crest SoftCall™ & JumboCall™ Momentary Call Cords

Crest SoftCall™ Momentary Call Cords

SoftCall provides a durable and affordable option for nurse call and light controls. The ergonomic pendant fits comfortably in the hand, and its silicone skin provides a reassuring, “soft” feel. The waterproof pendant can easily be wiped down for infection control. The high-contrast color makes the SoftCall pendant easy to find in bedding.


- Soft, silicone pendant
- Easy to clean
- Durable strain relief
- IP47 sealed & waterproof



LATEX-FREE

2
YEAR
WARRANTY

Part No.	System	Style	Color	Length	Price
SC44B-01 †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse	Blue	8'	54.39
SC44B-02 †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse	Blue	12'	56.70



Part No.	System	Style	Color	Length	Price
SC04B-01 †	Rauland Responder III+, IV, and 4000	Nurse	Blue	8'	48.32
SC04C-01 †	Rauland Responder III+, IV, and 4000	Nurse, 1 Light	Blue	8'	49.55
SC04E-01 †	Rauland Responder III+, IV, and 4000	Nurse, 2 Light	Blue	8'	50.56


Part No.	System	Style	Color	Length	Price
SC04C-02 †	Rauland Responder III+, IV, and 4000	Nurse, 1 Light	Blue	8'	49.59
SC04E-02 †	Rauland Responder III+, IV, and 4000	Nurse, 2 Light	Blue	8'	50.50

Crest JumboCall™ Momentary Call Cords

The JumboCall call cord is ideal for patients who have difficulty using their hands.

- Large, easy-to-grasp pendant
- For use in most systems using ¼" phone jacks

LATEX-FREE



Part No.	System	Style	Color	Length	From
5300G-__ †	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	Grey	10' or 12'	47.36


All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.




BreathCall® Air Activated Call Cords

BreathCall is for patients with little to no motor skills and is activated by simply breathing into a disposable straw and filter assembly. It has a flexible gooseneck with universal mounting clamp. The Mini BreathCall offers the same features of the standard with a more streamlined design to fit smaller areas. All BreathCalls come with a 9' cord, clamp and filter and straw kit. Standard BreathCalls have 36" black tubing from clamp to mouthpiece and Mini BreathCalls have 28" white tubing. **For single patient use.**

- Flexible conduit bends to almost any position to suit the user
- Sturdy clamp easily attaches to headboards
- Locking or non-locking switch closure available
- Also available with dual controls for nurse/TV or radio (11201-09D)



Part No.	System	Style	Color	Length	Price
11201-09	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single, Non-Locking	Black	9'	157.47
11201-09D	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	Black	9'	215.31
16001-09	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Mini, Single, Non-Locking	White	9'	122.26



Part No.	System	Style	Color	Length	Price
115576	Rauland Responder III+, IV, and 4000	Single, Non-Locking	Black	9'	181.86


All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V20. Revised 04/24/23.




Crest GeriCall™ Call Cords

The GeriCall Geriatric Call Cord is a round, flattened momentary call cord pad that works great for residents who have limited hand dexterity. This call cord utilizes a mechanical switch in the pad that can be activated by slight pressure from the hand, arm or body or can even be positioned under the resident's chin.

- Large pad for easy activation
- Includes security clip
- 3" diameter pad



Part No.	System	Style	Color	Length	From
118844C †	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single	Grey	8'	34.48
118844C-10 †	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single	Grey	10' or 15'	37.45
118844D-10 †	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double	Grey	10'	88.26




Part No.	System	Style	Color	Length	Price
118804C-10 †	Rauland Responder III+, IV, and 4000	Single	Grey	10'	51.97

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.

PressCall® Pneumatic Call Cords

PressCall is an air-activated call cord designed for patients who have difficulty activating standard call cords. PressCall has an oval-shaped, soft vinyl bulb design that fits easily into the hand.

- For use by patients with sufficient motor skills to press on the bulb
- Great for geriatric patients who cannot operate a regular call cord
- Easy to squeeze bulb style





LATEX-FREE

OXYGEN SAFE

1
YEAR LIMITED WARRANTY



1/4" Phone Plug
2-conductor





Part No.	System	Style	Color	Length	From
11372__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single, Non-Locking	White or Grey	6', 8', 10', 12' or 15'	39.91
1159__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	White or Grey	6', 8', or 12'	58.49
11373__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	White or Grey	8' or 10'	66.86

8-pin DIN plug




Part No.	System	Style	Color	Length	Price
113780†	Rauland Responder III+, IV, and 4000	Single, Non-Locking	Grey	8'	75.32

PadCall® Pneumatic Call Cords

PadCall is an air-activated call cord designed for patients who have difficulty activating standard call cords. PadCall has a large, circular, soft vinyl pad with gripping feet that prevent it from slipping off of any surface.

- For use by patients with sufficient motor skills to press or lean on the pad
- Great for geriatric patients who cannot operate a regular call cord
- Stay-put 3 1/2" wide pad





LATEX-FREE

OXYGEN SAFE

1
YEAR LIMITED WARRANTY



1/4" Phone Plug
2-conductor





Part No.	System	Style	Color	Length	From
1137__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single, Non-Locking	White or Grey	6', 8', 10', 12' or 15'	54.19
1159__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	White or Grey	6', 8', or 12'	77.98
1137__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	White or Grey	8' or 10'	98.61

8-pin DIN plug



Part No.	System	Style	Color	Length	Price
113781†	Rauland Responder III+, IV, and 4000	Single, Non-Locking	Grey	8'	86.39

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V20. Revised 04/24/23.



Crest TouchCall™ Call Cords

TouchCall provides reliable and versatile nurse call actuation for patients with limited dexterity. The TouchCall features a large, ultra-sensitive touch plate surface that can be easily pressed by a hand, foot, elbow, chin, etc. The tent design with clip and Velcro® attachment methods allow for multiple positioning possibilities to make it simple and convenient for any patient to use.



- Standard 10' cord length to suit a variety of applications
- Classic grey color for easy visibility

LATEX-FREE

OXYGEN SAFE



2
YEAR WARRANTY

1/4" Phone Plug
2-conductor



Part No.	System	Other Features	Style	Color	Length	Price
113673-4410†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Ultra-Sensitive Touch Plate	Single	Grey	10'	64.44

8-pin DIN plug



Part No.	System	Other Features	Style	Color	Length	Price
113673-0410†	Rauland Responder III+, IV, and 4000	Ultra-Sensitive Touch Plate	Single, Non-Locking	Grey	8'	71.98

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.




Crest Magnetic Break-A-Way Cable System

The Magnetic Break-a-Way Cable System saves both wall stations and pillow speakers from damage due to cable stress. Any cable pressure over 6 lb will break the magnetic connection and disconnect the pillow speaker without damage to the wall station. The magnets guide the connection, ensuring the pillow speaker is correctly and quickly connected - no more fumbling or broken pins.

- For use with Crest devices
- Eliminate downtime & cost of wall station repairs
- Save on pillow speaker replacements
- Cable length is approximately 24"
- Also works with call cords

LATEX-FREE



Part No.	System	Color	Length	Price
116075	Rauland Responder III+, IV, and 4000	White	2'	74.40




All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V20. Revised 04/24/23.


Cord Extensions & Station Savers

Cord Extensions provide an instant way to lengthen call cords for greater flexibility. Station Savers create a break-away point between cord and station to help prevent damage.

- Disconnection of cord triggers cord-out on stations with this feature
- Wired straight through for use with most 8-pin DIN phone plug systems
- Multiple lengths to fit any application



Part No.	System	Description	Color	Length	Price
114711-10	Rauland Responder III+, IV, and 4000	Cord Extension, right angle plug	White	10'	28.48
107397-10	Rauland Responder III+, IV, and 4000	Cord Extension	White	10'	39.95
107397	Rauland Responder III+, IV, and 4000	Station Saver	White	1'	25.03
114711-01	Rauland Responder III+, IV, and 4000	Station Saver, right angle plug	White	1'	25.27



Part No.	System	Description	Color	Length	Price
108090-01	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Station Saver, right angle plug	White	1'	24.07



Adapters

Type B Adapters convert your existing receptacle to the same receptacle and 1/4" phone receptacle. Include 9" extension cord and plug.

Type Y Adapters convert your existing receptacle to two 1/4" phone receptacles.

Type Y2 Adapters convert your existing receptacle to one female receptacle of the same style and one 1/4" phone receptacle.

Type D Adapters convert your existing receptacle to one 1/4" phone receptacle.

Part No.	System	Style	Description	Price
9515	Rauland Responder III+, IV, and 4000	Type B	8-pin DIN to 8-pin DIN and 1/4" 2C	81.32
CA04X00N2-01	Rauland Responder III+, IV, and 4000	Type Y	8-pin DIN to (2) 1/4" 2C	59.85
CA04S04N1-01	Rauland Responder III+, IV, and 4000	Type Y2	8-pin DIN to 8-pin DIN and 1/4" 2C	30.93
111541	Rauland Responder III+, IV, and 4000	Type D	8-pin DIN to 1/4" 2C	31.32

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.




B-32877P

B-32877

B-35509

Replacement Receptacles

Part No.	System	Description	Price
B-32877	Rauland Responder III+, IV, and 4000	8-pin DIN Receptacle, panel mount	7.95
B-32877P	Rauland Responder III+, IV, and 4000	8-pin DIN Receptacle, on 1-gang stainless steel plate	49.65
B-33509	Rauland Responder III+, IV, and 4000	8-pin DIN, Receptacle, PC mount	11.87



NCLED2CC


LED Corridor Lights

Crest LED replacement corridor lights for Rauland® Responder® IV NCLED2 and NCLED6. Can be mounted on a 1- or 2-gang box.

- Available in 2 or 6 lights: white, red, green, amber, yellow & blue
- Durable plastic lens

Part No.	System	Description	Price
NCLED2CC†	Rauland Responder IV	Corridor Light, 2-lamp, white and red	36.00
NCLED6CC†	Rauland Responder IV	Corridor Light, 6-lamp, white, green, amber, yellow, red and blue	42.18

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V20. Revised 04/24/23.




B-31993 Red

Dual Call Cord Patient Station

Crest replacement for Rauland® Responder® IV nurse call station NC2JACK features visual indicators for call-placed. Mounts on standard 1-gang back box, with a stainless steel faceplate.

- Calls can be reset at the station using the cancel pushbutton
- 14 VDC @ 100 mA



Part No.	System	Description	Price
NC2JACKC	Rauland Responder IV	Patient station, dual ¼" jacks	105.01
B-31993	N/A	Replacement Bulb, red	0.18



Auxiliary Alarm 4 Input Station

Crest replacement for Rauland® Responder® IV and Responder® 4000 nurse call station NC4JACK includes four ¼" 2-conductor jacks. Mounts on standard 1-gang back box, with a stainless steel faceplate.

Part No.	System	Description	Price
NC4JACKC	Rauland Responder IV/4000	Patient station, ¼" jacks	63.04



NCPCS1C


NCPBPHC

Pushbutton & Waterproof Pullcord Stations


Crest replacement Pushbutton and Waterproof Pullcord Stations for Rauland® Responder® IV system. Mounts on a 1-gang box, with a stainless steel or waterproof plastic faceplate. Call-placed light included.

Part No.	System	Description	Price
NCPCS1C†	Rauland Responder IV	Pullcord Bath Station, replacement for NCPCS1, waterproof, with white CleanCord	171.96
NCPBPHC	Rauland Responder IV	Pushbutton Station, replacement for NCPBPH	150.05
B-31993	N/A	Replacement Bulb, red	0.18


All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.



R4K11VC



QUALITY REPLACEMENT




B-31993 Red

Patient Stations


Crest replacement for Rauland® Responder 4000 Visual Call Cord Stations R4K11V and R4K21V. Single and dual stations feature a red call-placed indicator light and come standard with 2-gang stainless steel plate that mounts on a 1-gang electrical box.

- Calls can be reset at the station using the cancel pushbutton
- Accepts 8-pin DIN call cords
- 15.5 VDC @ 0.018A


Part No.	System	Description	Price
R4K11VC	Rauland Responder 4000	Patient Station, single	196.68
R4K21VC	Rauland Responder 4000	Patient Station, dual	214.19
B-31993	N/A	Replacement Bulb, red	0.18



R4K22AC



QUALITY REPLACEMENT




B-31993 Red
E-11657 Green

Patient Stations


Crest replacements for Rauland® Responder 4000 patient stations for use with 8-pin DIN call cords. Single and dual stations feature visual indicators for call placed/cord out (red) and communication/intercom on (green).

- Calls can be reset at the station using the cancel pushbutton
- 14 VDC @ 100 mA
- 4.5" x 4.562" 2-gang stainless steel plate


Part No.	System	Description	Price
R4K12AC †	Rauland Responder 4000	Patient Station, single	135.50
R4K22AC †	Rauland Responder 4000	Patient Station, dual	163.72
B-31993	N/A	Replacement Bulb, red	0.18
E-11657	N/A	Replacement Bulb, green	



R4KPC10C



QUALITY REPLACEMENT




B-31993 Red

Pullcord Station


Non-waterproof Crest replacement for Rauland® Responder 4000 nurse call station. Mounts on a 2-gang box stainless steel faceplate.

Part No.	System	Description	Price
R4KPC10C †	Rauland Responder 4000	Pullcord Station, with white CleanCord	126.19
B-31993	N/A	Replacement Bulb, red	0.18


All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 02/03/21.




BS800C




BS810C




BS820C



BS300C



QUALITY REPLACEMENT



B-31993 Red
E-11657 Green

Patient Station Replacements for Rauland® Responder® III+

Crest replacement stations for Rauland Responder III+ nurse call systems feature visual indicators for call placed and station power. Mounts on a 1-gang box, with a stainless steel faceplate.

- Calls can be canceled at the stations using reset button
- Available in ¼" or 8-pin DIN styles


Part No.	System	Description	Price
BS800C †	Rauland Responder III+	Bedside Station, dual ¼" jacks	57.57
BS810C †	Rauland Responder III+	Patient Station, single, 8-pin DIN jack	82.70
BS820C †	Rauland Responder III+	Bedside Station, dual 8-pin DIN jacks	162.87
BS300C †	Rauland Responder III+	Bedside Station with Speaker, single, ¼" jack	144.13
B-31993	N/A	Replacement Bulb, red	0.18
E-11657	N/A	Replacement Bulb, green	

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V20. Revised 04/24/23.



Facility & Clinical Management Solutions

Crest is your *First Call* for all of your *Critical Needs*



PHONE: 1-800-328-8908 | FAX: 1-800-369-9207 | WWW.CRESTHEALTHCARE.COM

FULL CATALOG AVAILABLE ONLINE

www.CrestHealthcare.com

- Over 3,000 products available
- Order all your facility products
- Bookmark, print or email pages



Ready to Place An Order?

Give us a call or visit us online to get started today.

1-800-328-8908

www.CrestHealthcare.com