



**Crest Healthcare**

---

# **TekTone® Tek-CARE®300III Parts List**

Your guide to TekTone® Tek-CARE®300III system parts available from Crest that are durable, economical and designed with infection control in mind.

**CREST**  
HEALTHCARE®

© 2023 Crest Healthcare





## Enhancing the environment of care for those in need and those who care

Crest Healthcare is a leading provider of a wide range of healthcare products and technical services. Historically specialized in high-quality nurse call replacement parts, and today producing nurse call system components for nearly every current system, Crest has been a trusted source for innovative product development and production since 1967.

With customer satisfaction as our driving goal, Crest products provide quality, value, innovation, variety and durability. Working together to guarantee a positive experience, our in-house team includes customer service, technical support, engineering, manufacturing and quality assurance experts.

### Engineering

Our engineering team is continually creating and testing new communication products. With the use of destructive and non-destructive in-house testing and prototype evaluations, Crest ensures every product provides durability and a long service life. Custom services are available to OEMs and other commercial customers.

### Manufacturing

We build and assemble all the Crest products in-house and work closely with our manufacturing partners to ensure our same high quality standards are met.

### Quality Assurance

We conduct testing on 100% of finished products we manufacture prior to shipment. Crest products meet or exceed demanding safety standards up to and including UL 1069 Recognition.

### Quality Certifications

Crest strives to offer you high quality products to support your facility needs. Crest Healthcare is ISO 13485:2016 certified provider of Pillow Speakers, Nurse Call Stations, Nurse Call Cords, and Bed Control Equipment. Our certification represents our commitment to on-going quality improvements and a dedication to higher levels of performance, service and quality products. We received our first certification in 1999 (ISO 9001) and Crest continues to improve quality in all aspects of designing, manufacturing and customer service.

# TekTone® Tek-CARE®300III

## Nurse Call Components

### Pillow Speakers

Elite Digital™ .....	4
Elite Infrared™ .....	5
Replacement Cable Assemblies .....	6
Remote Call-IR™ .....	7
Slimline™ .....	8

### Call Cords

#### Momentary Call Cords

SoftCall™ .....	9
EconoCall™ .....	10
CleanCall™ .....	11
DuraCall™ .....	12
DuraColor™ .....	13
DuraLight™ .....	13
DuraGlo™ .....	13

#### Specialty Call Cords

GeriCall® .....	14
BreathCall™ .....	15
PressCall® .....	16
PadCall® .....	17
TouchCall™ .....	18

### Additional Components

Magnetic Break-A-Way Cable System .....	19
Call Cord Extenders .....	19
Station Savers .....	20
Dummy Plugs .....	20
Adapters .....	21

### Stations

Visual Signals .....	22 - 23
Wall Stations .....	24 - 29

## EliteDigital™

### Crest Elite Digital™ Pillow Speakers

Designed for use with digital TVs, Autosense™ Technology automatically identifies the brand of digital healthcare TV and sends the correct signal. They can be moved to any room, regardless of healthcare television brand.

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille

#### Specifications

Height:	9"
Width:	4 3/5"
Thickness:	1 3/4"
Cord Length:	8'
Case Material:	ABS, UL94V-0
Cord Connector:	Plug-and-Play (solderless)
Switch Rating:	30 VDC @ 100 mA
Warranty:	2-year warranty against defects in materials and workmanship from date of purchase

LATEX-FREE

OXYGEN SAFE



QUANTITY DISCOUNTS  
CALL FOR PRICING

2  
YEAR WARRANTY



1/4" Phone Plug  
3-conductor



Part No.	System	Nurse Call	TV Control	Other Features	Price
EN0301EEAAXXXAE1†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter		Up/Down		112.69
EN0301EXAAXXXAE1†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter		10-key		112.63

8-pin DIN  
Plug



Part No.	System	Nurse Call	TV Control	Other Features	Price
XN1601EBA500A10N1†	Tek-CARE300III	Nurse	Up/Down		135.39
XN1601EEA500A10N1†	Tek-CARE300III	Nurse	Up/Down	2 lights	129.25
XN1601XBA500A10N1†	Tek-CARE300III	Nurse	10-key		119.95

RJ45 8-pin  
Plug



Part No.	System	Nurse Call	TV Control	Other Features	Price
EN1301EXABA200AE1†	Tek-CARE300III	Nurse	10-key		119.66
EN1302EXABA200ZG1†	Tek-CARE300III	Nurse	10-key	Headphone	142.92





# EliteInfrared™

## Crest Elite Infrared™ Pillow Speakers

The Crest Elite Infrared Pillow Speaker is the only pillow speaker that combines nurse call and TV with external DVD or Set Top Box (STB) controls. It combines the ease of a universal remote with traditional pillow speaker functionality such as nurse call and environmental controls. It is also useful where there is no cross-room wiring or with TVs that do not have a pillow speaker port. Uses two AAA batteries (included).


- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille

### Specifications

<b>Height:</b>	9"
<b>Width:</b>	4 ⅜"
<b>Thickness:</b>	1 ¾"
<b>Cord Length:</b>	8'
<b>Case Material:</b>	ABS, UL94V-0
<b>Cord Connector:</b>	Plug-and-Play (solderless)
<b>Power Requirements:</b>	2 AAA alkaline batteries required (included) (Infrared only)
<b>Switch Rating:</b>	30 VDC @ 100 mA
<b>Warranty:</b>	2-year warranty against defects in materials and workmanship from date of purchase




¼" Phone Plug  
3-conductor



Part No.	System	Nurse Call	TV Control	Other Features	Price
ET0301EXAAXXXZE1†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter		10-key		<b>133.16</b>
ET0301EXAAXXXZG1†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter		10-key	Headphone	<b>150.20</b>

¼" Phone Plug  
2-conductor



Part No.	System	Nurse Call	TV Control	Other Features	Price
ET4401EXABA600ZE1†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse	10-key		<b>146.80</b>
ET4401EDABA600ZE1†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse	10-key	DVD	<b>164.46</b>

8-pin DIN  
Plug



Part No.	System	Nurse Call	TV Control	Other Features	Price
XT1601EBA500Z10N1†	Tek-CARE300III	Nurse	Up/Down		<b>142.57</b>
XT1601XEA500Z10N1†	Tek-CARE300III	Nurse	10-key	2 lights	<b>148.46</b>

RJ45 8-pin  
Plug



Part No.	System	Nurse Call	TV Control	Other Features	Price
ET1301EXABA400ZE1†	Tek-CARE300III		10-key		<b>148.59</b>
ET1301EEABA200ZE1†	Tek-CARE300III	Nurse	Up/Down		<b>141.00</b>

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.

# Replacement Cable Assemblies



XN4401W-08



XN1601W-08



XN1301W-08

## Replacement Cable Assemblies for Elite Pillow Speakers

Replacement pillow speaker cable assemblies for Elite Digital™ and Elite Infrared™ are easy to replace using "plug-and-play" connectors. Standard white cables are 8' with molded-in strain relief. To verify the correct cable for your pillow speaker or for replacement assistance, please contact Crest Technical Support.

1/4" Phone Plug  
2-conductor



Part No.	System	Color	Length	Price
XN4401W-08	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	White	8'	<b>25.95</b>

1/4" Phone Plug  
3-conductor



Part No.	System	Color	Length	Price
XN0301W-08	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	White	8'	<b>31.06</b>

8-pin DIN  
Plug



Part No.	System	Color	Length	Price
XN1601W-08	Tek-CARE300III	White	8'	<b>34.97</b>

RJ45 8-pin  
Plug



Part No.	System	Color	Length	Price
XN1301W-08	Tek-CARE300III	White	8'	<b>34.91</b>

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.





Programming instructions and code list can be found at [www.cresthealthcare.com](http://www.cresthealthcare.com) under the "Resources" tab

## Crest Remote Call-IR™ TV Control & Nurse Call

Replace your call cord and remote control with one simple device. The Crest Remote Call-IR controls the TV using IR signals like a universal remote while still providing the safety of reliable wired nurse call.

Remote Call-IR plugs into your nurse call station to ensure it stays within the room and is durable enough to withstand drops and daily use. The 10-key overlay style is great for emergency rooms or private patient rooms, where the room is not wired for TV sound or the sound control of a pillow speaker is not needed. The large button overlay is perfect for kids, geriatric patients or residents with limited dexterity. Features simplified control of TV power, channel and volume with big, easy-to-see buttons. Both styles feature an on-board library of codes that can be used to program it to work with almost any TV. Uses two AA batteries (not included) that last as long as in typical TV remotes.

- Robust design, easy to disinfect
- No batteries needed for nurse call operation
- Distinctive nurse call button with Braille
- Bed A & B functionality
- Fits Crest Pillow Speaker Keeper

¼" Phone Plug  
2-conductor



## Specifications

<b>Height:</b>	6 ½"
<b>Width:</b>	Tapers, 2 ⅝" - 2 ¼"
<b>Thickness:</b>	1 ⅜" at thickest point
<b>Weight:</b>	13 ⅞ oz (with batteries)
<b>Cord Length:</b>	10'
<b>Case Material:</b>	ABS, Flame rated (UL94V-0)
<b>Cord &amp; Strain Relief:</b>	Flexible molded PVC
<b>Cord Connector:</b>	Plug-and-Play (solderless)
<b>Plug Style:</b>	2-conductor or 3-conductor ¼" phone plug (call Crest for available adapters)
<b>Power Requirements:</b>	Two AA alkaline batteries (sold separately)
<b>Switch Rating:</b>	100 mA @ 28 VDC
<b>Warranty:</b>	2-year warranty against defects in materials and workmanship from date of purchase

LATEX-FREE

OXYGEN  
SAFE

QUANTITY  
DISCOUNTS  
CALL FOR PRICING

2  
YEAR  
WARRANTY



114576

Part No.	System	Nurse Call	Style	Length	Price
XP441E0B600XXX2†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse	Large Button	10'	83.58
XP441X0B600XXX2†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse	10-key	10'	89.95

## Replacement Cable for Remote Call IR

Part No.	System	Length	Price
XP4401W-10	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	10'	22.95

¼" Phone Plug  
3-conductor



Part No.	System	Nurse Call	Style	Length	Price
XP031E0AXXXXXX2	Any station with ¼" jack or Type D, B, Y, or Y2 adapter		Large Button, Headphone	10'	86.25

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.





# Slimline™

## Crest Slimline™ Pillow Speakers

Slimline Pillow Speakers are designed for each hospital's specific needs and fit multiple nurse call systems. These units feature a slim "easy to hold" design, with many options available. Slimline pillow speakers are priced far below original equipment cost and are designed to make repairs quick and easy.

- Slim, easy to hold design
- Comes standard with 8' cord length

RJ45 8-pin Plug



### Specifications

Height:	7.5"
Width:	2.6" (tapering to 1.9")
Thickness:	1.7"
Weight:	18.2 oz
Cord Length:	8'
Color:	White
Case Material:	UL 94V-0 Flame Rated, Latex Free ABS
Speaker Material:	Moisture-resistant mylar cone
Electrical Rating:	24 VDC at 100 mA
Standards:	Tested to UL 1069 O2 Safe Testing with oxygen enriched environment marking RoHS Compliant FDA Listed
Warranty:	2-year warranty against defects in materials and workmanship from date of purchase



Part No.	System	Nurse Call	TV Control	Other Features	Price
PSS13010801†	Tek-CARE300III	No	Yes		87.67
PSS13020801†	Tek-CARE300III	Yes	Yes		82.57
PSS13050801†	Tek-CARE300III	Yes	Yes	2 Lights	89.30

### Replacement Cable Assemblies for Crest Slimline™

Replacement cable assemblies for Slimline are easy to replace. Standard white cables are 8' with molded-in strain relief. To verify the correct cable for your pillow speaker or for replacement assistance, please call Crest Technical Support.

#### Slimline Cable Assemblies

Part No.	Style	System	Color	Length	Price
PS1301W-08	RJ45 8-pin	Tek-CARE300III	White	8'	18.24

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.





LATEX-FREE


2  
YEAR  
WARRANTY

### Crest SoftCall™ Waterproof Nurse Call & Light Control

SoftCall provides a durable and affordable option for nurse call and controls. The ergonomic pendant fits comfortably in the hand, and its silicone skin provides a reassuring, "soft" feel. The waterproof pendant can easily be wiped down for infection control. The high-contrast color makes the SoftCall pendant easy to find in bedding.

- Soft silicone pendant
- Easy to clean
- Durable strain relief
- IP47 sealed & waterproof
- Standard 8' white cable

¼" Phone Plug  
2-conductor



Part No.	System	Nurse Call	Other Features	Style	Color	Length	Price
SC44B-01	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse		Single	Blue	8'	<b>43.55</b>
SC44B-02	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse		Single	Blue	12'	<b>45.42</b>

RJ45 8-pin  
Plug



Part No.	System	Nurse Call	Other Features	Style	Color	Length	Price
SC13B-01	Tek-CARE300III	Nurse		Single	Blue	8'	<b>57.56</b>
SC13C-01	Tek-CARE300III	Nurse	1 Light	Single	Blue	8'	<b>58.85</b>
SC13E-01	Tek-CARE300III	Nurse	2 Lights	Single	Blue	8'	<b>59.71</b>

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.

# Crest EconoCall™ Momentary Call Cords



## Crest EconoCall™ Call Cords

The EconoCall delivers a clean, basic design with reliable function at a price any facility can afford.

- Our most economically priced call cord
- Repairable button for fast fixes
- Plug styles available to fit most systems
- RoHS 2 Compliant

Straight  
¼" Phone Plug  
2-conductor



Part No.	System	Style	Color	Length	From
117700-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	White	7' or 12'	18.09

¼" Phone Plug  
2-conductor



Part No.	System	Style	Color	Length	From
7700W-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	White	7', 10', 12', 15' or 20'	13.77
7700B-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	Beige	7', 10', or 12'	16.69
7700G-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Single	Grey	7', 10' or 12'	18.83
7700DW-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Double	White	7', 10', 12', 15' or 20'	24.71
7700DB-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Double	Beige	10' or 12'	31.54
7700DG-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Double	Grey	10' or 12'	31.54

8-pin DIN  
Plug



Part No.	System	Style	Color	Length	Price
117716†	Tek-CARE300III	Single	White	7'	21.91
117716-10†	Tek-CARE300III	Single	White	10'	29.26

RJ45 8-pin  
Plug



Part No.	System	Style	Color	Length	From
7710W-__†	Tek-CARE300III	Single	White	7', 10', 12', 15' or 20'	20.54
7710DW-12†	Tek-CARE300III	Double	White	12'	47.64

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.





## Crest CleanCall™ Call Cords

Maintaining a clean environment can be difficult, but CleanCall simplifies the job. The completely sealed membrane style button can be easily wiped down to prevent the spread of germs, making it a great option for infection control in any facility.

- Smooth membrane is easy to sanitize
- Easy grasp winged design fits comfortably in the hand
- Highly visible, high-contrast button with nurse call symbol
- RoHS 2 Compliant


LATEX-FREE

OXYGEN  
SAFE

1

YEAR  
WARRANTY

1/4" Phone Plug  
2-conductor



Part No.	System	Style	Color	Length	From
8900W-__†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single	White	8', 10', 12' or 15'	<b>27.26</b>
8900DW-10†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double	White	10'	<b>61.97</b>

8-pin DIN  
Plug



Part No.	System	Style	Color	Length	Price
8916W-8†	Tek-CARE300III	Single	White	8'	<b>30.59</b>

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.

# Crest DuraCall™ Momentary Call Cords



## Crest DuraCall™ Call Cords

When durability and feature versatility are a must, look no further than the DuraCall Call Cord line. With its rugged design, DuraCall has survived drops, submersion in water and even being driven over by a car.

- Easy-grip winged pendant
- Tactile "click" raised button
- Easy to clean
- Wide variety of plug styles & call cord features available
- Short 6" length available for safety in anti-ligature environments



Durable Waterproof Design

1/4" Phone Plug  
2-conductor



Part No.	System	Style	Color	Length	From
9900W-.5†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Anti-Ligature	White	6"	29.38
9900W-__†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single	White	7', 10', 12', 15' or 20'	16.05
9900DW-__†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double	White	7', 10', 12', 15' or 20'	25.50
9900CLW-16†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single, Coiled	White	16'	29.08

8-pin DIN  
Plug

Part No.	System	Style	Color	Length	From
9916W-__†	Tek-CARE300III	Single	White	7' or 10'	23.91

RJ45 8-pin  
Plug

Part No.	System	Style	Color	Length	From
9910W-__†	Tek-CARE300III	Single	White	7' or 10'	21.09
9910DW-7†	Tek-CARE300III	Double	White	10'	28.67

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.




## Crest DuraColor™ Call Cords

Crest's neon green pendant in the DuraCall Call Cord line is the latest option for nurse call visibility and ease of use. Bright and vibrant, the pendant is ideal for those with poor eyesight or for those who may have difficulty distinguishing it from bedding.

- High contrast neon color for easy identification
- Prominent color assists patients & staff
- White cable



¼" Phone Plug  
2-conductor



LATEX-FREE



2  
YEAR  
WARRANTY

Part No.	System	Other Features	Style	Color	Length	From
9900N-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	High-Contrast Color Pendant	Single	Green	7' or 10'	21.68


## Crest DuraLight™ Call Cords

DuraLight Call cords for ¼" phone plug systems have LEDs that momentarily light up when the button is depressed (may not work with every system).

- Bright red LEDs light up the pendant
- Requires no batteries
- FDA Listed



¼" Phone Plug  
2-conductor



LATEX-FREE



2  
YEAR  
WARRANTY

Part No.	System	Other Features	Style	Color	Length	From
9900R-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Red LEDs	Single	White	7' or 10'	32.52


## Crest DuraGlo™ Call Cords

DuraGlo is the only call cord in the industry with glowing bands to help your residents locate their call cord in dark rooms! DuraGlo can glow up to four hours in total darkness. Built on the same solid platform as our popular DuraCall call cord. Works with any nurse call system.

- Easily locate in dark rooms
- Recharges with light exposure
- Requires no power or batteries



¼" Phone Plug  
2-conductor



LATEX-FREE



2  
YEAR  
WARRANTY

Part No.	System	Other Features	Style	Color	Length	From
9900L-__†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Glow-in-the-Dark Pendant	Single	White	7' or 10'	33.29

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.



Crest GeriCall™ Call Cords

The GeriCall Geriatric Call Cord is a round, flattened momentary call cord pad that works great for residents who have limited hand dexterity. This call cord utilizes a mechanical switch in the pad that can be activated by slight pressure from the hand, arm or body, or even can be positioned under the resident's chin.

- Large pad for easy activation
- Includes security clip
- 3" diameter pad




LATEX-FREE

OXYGEN SAFE


2 YEAR WARRANTY

1/4" Phone Plug  
2-conductor




Part No.	System	Style	Color	Length	From
118844C†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single	Grey	8'	32.92
118844C-__†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single	Grey	10' or 15'	35.10
118844D-10†	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double	Grey	10'	88.26

8-pin DIN  
Plug



Part No.	System	Style	Color	Length	Price
118816C-10†	Tek-CARE300III	Single	Grey	10'	50.34

RJ45 8-pin  
Plug



Part No.	System	Style	Color	Length	Price
118813C-10†	Tek-CARE300III	Single	Grey	10'	55.43

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.



11201-09

## BreathCall® Air Activated Call Cords

BreathCall is for patients with little to no motor skills and is activated by simply breathing into the disposable straw and filter assembly. It has a flexible gooseneck with universal mounting clamp. The Mini BreathCall offers the same features of the standard with a more streamlined design to fit smaller areas. All BreathCalls come with a 9' cord, clamp and filter and straw kit. Standard BreathCalls have 36" black tubing from clamp to mouthpiece, and Mini BreathCalls have 28" white tubing. For single patient use.

- Flexible conduit bends to almost any position
- Sturdy clamp easily attaches to headboards
- Locking or non-locking switch closure available
- Dual controls for nurse, TV or radio available (11201-09D)



LATEX-FREE

OXYGEN  
SAFE

1

YEAR LIMITED  
WARRANTY

¼" Phone Plug  
2-conductor



11001-WE8-09

Part No.	System	Other Features	Style	Color	Length	Price
11201-09	Any station with ¼" jack or Type D, B, Y, or Y2 adapter		Single, Non-Locking	Black	9'	<b>167.30</b>
11201-09D	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Nurse, TV or Radio controls	Double, Non-Locking	Black	9'	<b>220.86</b>
16001-09	Any station with ¼" jack or T type D, B, Y, or Y2 adapter		Mini, Single, Non-Locking	White	9'	<b>133.67</b>

RJ45 8-pin  
Plug



Part No.	System	Other Features	Style	Color	Length	Price
11001-WE8-09	Tek-CARE300III		Single, Non-Locking	Black	9'	<b>203.52</b>

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.



PressCall® Pneumatic Call Cords

PressCall is a pneumatic air-activated call cord designed for patients who have difficulty activating standard call cords. The PressCall has an oval shaped, soft vinyl bulb design that fits easily in the hand.

- For patients with sufficient motor skills too press bulb
- Great for patients unable to operate regular call cords
- Easy to squeeze bulb style












1/4" Phone Plug  
2-conductor






Part No.	System	Style	Color	Length	From
11372_	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single, Non-Locking	White or Grey	6', 8', 10', 12' or 15'	42.65
1159__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	White or Grey	6', 8', or 12'	60.85
11373_	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	White or Grey	8' or 10'	67.58


8-pin DIN  
Plug






Part No.	System	Style	Color	Length	Price
115578	Tek-CARE300III	Single, Non-Locking	Grey	8'	80.71

RJ45 8-pin  
Plug





Part No.	System	Style	Color	Length	Price
113782	Tek-CARE300III	Single, Non-Locking	Grey	8'	82.85

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.





## PadCall® Pneumatic Call Cords

PadCall is a pneumatic air-activated call cord designed for patients who have difficulty activating standard call cords. They have a large circular, soft vinyl pad with gripping feet that prevent it from slipping off any surface.

- For patients with sufficient motor skills too press or lean on pad
- Great for patients unable to operate regular call cords
- Stay-put 3 1/2" wide pad




LATEX-FREE



1  
YEAR LIMITED  
WARRANTY

1/4" Phone Plug  
2-conductor



Part No.	System	Style	Color	Length	From
1137__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Single, Non-Locking	White or Grey	6', 8', 10', 12' or 15'	61.76
1159__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	White or Grey	6', 8', or 12'	91.28
1137__	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	Double, Non-Locking	White or Grey	8' or 10'	94.74

8-pin DIN  
Plug



Part No.	System	Style	Color	Length	Price
115579	Tek-CARE300III	Single, Non-Locking	Grey	8'	90.46

RJ45 8-pin  
Plug



Part No.	System	Style	Color	Length	Price
113783	Tek-CARE300III	Single, Non-Locking	Grey	8'	95.79

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.



## Crest TouchCall™ Call Cords

TouchCall provides reliable and versatile nurse call actuation for patients with limited dexterity. The TouchCall features a large, ultra-sensitive touch plate surface that can be easily pressed by hand, foot, elbow, chin, etc. The tent design with clip and Velcro® attachment methods allow for multiple positioning possibilities to make it simple and convenient for any patient to use.


- Standard 10' cord length to suit a variety of applications
- Classic grey color for easy visibility

LATEX-FREE

OXYGEN  
SAFE

2  
YEAR  
WARRANTY

¼" Phone Plug  
2-conductor



Part No.	System	Other Features	Style	Color	Length	Price
113673-4410†	Any station with ¼" jack or Type D, B, Y, or Y2 adapter	Ultra-Sensitive Touch Plate	Single	Grey	10'	57.27

8-pin DIN  
Plug



Part No.	System	Other Features	Style	Color	Length	Price
113673-1610†	Tek-CARE300III	Ultra-Sensitive Touch Plate	Single	Grey	10'	59.60

RJ45 8-pin  
Plug



Part No.	System	Other Features	Style	Color	Length	Price
113697-1310†	Tek-CARE300III	Ultra-Sensitive Touch Plate	Single	Grey	10'	61.91

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. † Item is FDA listed. V12. Revised 11/28/23.



## Magnetic Break-A-Way Cable System

The Magnetic Break-a-Way Cable System saves both wall stations and pillow speakers from damage due to cable stress. Any cable pressure over 6 lb will break the magnetic connection and disconnect the pillow speaker without damage to the wall station. The magnets guide the connection, ensuring the pillow speaker is correctly and quickly connected - no more fumbling or broken pins.

- Eliminate downtime & cost of wall station repairs
- Save on pillow speaker replacements
- Cable length is approximately 2'
- Also works with call cords

LATEX-FREE



8-pin DIN  
Plug



Part No.	System	Color	Length	Price
116075	Tek-CARE300III	White	2'	54.11

RJ45 8-pin  
Plug



Part No.	System	Color	Length	Price
116099	Tek-CARE300III	White	2'	59.94

## Call Cord Extenders

Call Cord Extenders allow you to lengthen any 2-conductor, 1/4" phone plug nurse call cord.

- Disconnection of cord triggers cord-out on stations with this feature
- Multiple lengths to fit any application

1/4" Phone Plug  
2-conductor



108090


Part No.	System	Color	Length	Price
108090-01	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	White	1'	38.50
108090	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	White	8'	24.03
108090-12	Any station with 1/4" jack or Type D, B, Y, or Y2 adapter	White	12'	24.55

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.

# Station Savers & Dummy Plugs



114681-10



114681-01

## Station Savers

Station Savers create a break-away point between cord and station to help prevent damage. Disconnection of cord triggers cord-out on stations with this feature.

Part No.	System	Description	Length	Price
114681-01	Tek-CARE300III	RJ45 8-pin	1'	22.81
114681-10	Tek-CARE300III	RJ45 8-pin	10'	28.53



113703



7710W-D

## Dummy Plugs

Dummy plugs provide proper termination to eliminate cord out calls from being placed when a call cord, pillow speaker or bed is disconnected from the wall station. Select the plug end type to fit your needs.

Part No.	System	Description	Price
113703	Tek-CARE300III	8-pin DIN	11.27
7710W-D	Tek-CARE300III	RJ45 8-pin, 5" connecting cord	14.83

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.



## Crest Nurse Call Adapters

Adapters accommodate multiple nurse call devices from a single wall receptacle or convert your existing station to an alternative plug style. Allows patient monitoring equipment to be used simultaneously with a standard call cord or pillow speaker.

**Type B Adapters** convert your existing receptacle to the same receptacle and one ¼" phone receptacle.

**Type D Adapters** convert your existing receptacle to one ¼" phone receptacle.

**Type Y Adapters** convert your existing receptacle to two ¼" phone receptacles.

**Type Y2 Adapters** convert your existing receptacle to one female receptacle of the same style and one ¼" phone receptacle.

**Type Y3 Adapters** convert your existing receptacle to two of your existing receptacles.



CA16S16N1-01



CA16D16N1-01



CA16X00N1-01

Part No.	System	Style	Description	Price
CA16S16N1-01	Tek-CARE300III	Type Y2	8-pin DIN to 8-pin DIN and ¼" 2-conductor	57.11
CA16D16N1-01	Tek-CARE300III	Type Y3	8-pin DIN to (2) 8-pin DIN	37.67
CA16X00N1-01	Tek-CARE300III	Type D	8-pin DIN to ¼" 2-conductor	33.66



CA13X00N2-01



9514



CA13S13N1-01



108058



115917

Part No.	System	Style	Description	Price
CA13X00N2-01	Tek-CARE300III	Type Y	RJ45 8-pin to (2) ¼" 2-conductor	57.11
CA13S13N1-01	Tek-CARE300III	Type Y2	RJ45 8-pin to RJ45 8-pin and ¼" 2-conductor	66.18
108058	Tek-CARE300III	Type D	RJ45 8-pin to ¼" 2-conductor	31.83
115917	Tek-CARE300III	Type D	RJ45 8-pin to ¼" 3-conductor	29.20
9514	Tek-CARE300III	Type B	RJ45 8-pin to RJ45 8-pin and ¼" 2-conductor	96.97

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.




4-Lamp LED Corridor Dome Light – OEM

The TekTone 4-lamp LED Corridor Dome Light provides visual indication that a call was made. Amber, green, red and white colored LEDs indicate patient call priority and staff presence. When used with any room station or IR310E TekTone Staff Station, the LI380LED indicates calls originated from that specific room. Fits standard 1-gang or 2-gang box with 1-gang or 2-gang ring.

- For use with Tek-CARE300III System
- Long lasting LEDs, no bulb replacement required
- Plug-in connector included
- Attractive flame-retardant plastic panel
- UL® 1069 Listed

Part No.	System	Description	Price
LI380LED	Tek-CARE300III	LED Corridor Pilot Light, 4-lamp, 24 VDC	85.17



Addressable Dome/Zone Light – OEM

TekTone's LI386LED Dome/Zone Light with built in multipurpose station provides visual indication of calls originating from bath, emergency and code call switches when used as a dome light. The LI386LED indicates calls originating from patient and staff stations when used as a zone light. The LI386LED also indicates all calls from the zones to which those stations are assigned and provides communication to and from the central equipment, including the nurse master station and monitor. No controls are provided.

- For use with Tek-CARE300III System
- Long lasting LEDs
- No bulb replacement required
- Simplified multiplex wiring
- UL® 1069 Listed

Part No.	System	Description	Price
LI386LED	Tek-CARE300III	LED Corridor Zone Light, 24 VDC	283.90

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.



LI381




### Corridor Dome Lights – OEM

TekTone original equipment Corridor Dome Lights provide visual indication of any type call and fit standard 2-gang, plaster ring. Features translucent white plastic cover which mounts on a standard mounting plate.

- For use with Tek-CARE300III System
- Designed for easy bulb replacement
- UL®1069 Listed

Part No.	System	Description	Price
LI381	Tek-CARE300III	Corridor Dome Light, 1-lamp, 28 VDC, 40 mA	26.15
LI382	Tek-CARE300III	Corridor Dome Light, 2-lamp, 28 VDC, 40 mA	40.08
RP092		Replacement Lens for Corridor Dome Light	7.97



### LED Corridor Zone Light – OEM

The LI382LED Corridor Zone Light provides visual indication of calls when used with a nurse call system, emergency call system, or similar application. Indicates an emergency call by the rapid flashing of the two red LED bulbs associated with the zone or area from which an emergency call has been placed. Normal calls are indicated by a steady illumination of the two white LEDs. If an emergency or normal call are placed in the same zone at the same time, the red LEDs connected to the emergency station will flash rapidly, while the white LEDs associated with the normal call will maintain a steady illumination.

- For use with Tek-CARE300III System
- Long lasting LEDs & no bulb replacement needed
- Plug-in connector supplied
- Attractive flame-retardant plastic panel
- Standard mounting plate for 1-gang or 2-gang mounting options
- UL®1069 Listed

Part No.	System	Description	Price
LI382LED	Tek-CARE300III	LED Corridor Zone Light, 24 VDC	66.80

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.



NC415G3  
Shown with optional TA415L



Master Station with LCD Touchscreen Display - OEM

This master station console option for systems using Gen2 central equipment modules provides for push to talk audio control for audio capable systems, but if desired, add the optional handset kit for better privacy. The user interface layout allows staff to view calls on the screen across devices on the same system. Up to three incoming calls are displayed by priority and time and are also color coded for easy recognition. Routine calls are annunciated with a steady call LED and slow pulsing tone with station number displayed on the touchscreen. Features a color LCD display with touch screen, with controls for tone volume, page view buttons, clock, and call listing. All calls must be canceled from the point of origin unless it is a routine call from an audio capable station.

- For use with Tek-CARE300III System
- Touchscreen to simply navigate calls or menus
- Use for both audio/visual & tone/visual applications
- Calls can be selected out of order
- RoHS compliant

Part No.	System	Description	Price
NC415G3	Tek-CARE300III	Master Station with LCD Touchscreen Display	1460.03
TA415L	Tek-CARE300III	Optional Handset Kit, includes cradle, handset and cord	200.00
IH415LD	Tek-CARE300III	Optional Desk Mount Housing	25.05
IH415W	Tek-CARE300III	Optional Wall Mount Bracket	33.40

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.





IR311E



IR319B



IR320B




## Patient Stations – OEM

The IR319B Single Patient Station and the IR320B Dual Patient Station from TekTone have two-way communication between patient and nurse, and provides for programmable privacy, paging and patient priority from the master station. If a call cord or pillow speaker is accidentally removed from its receptacle, the nurse master station registers a "cord out" call. Dual Patient Station additionally features independently operating receptacle and a second call-placed indicator.

- For use with Tek-CARE300III System
- Select 8-pin or RJ45 receptacles
- LED light indicate call-placed & intercom in use
- UL® 1069 Listed

Part No.	System	Description	Price
IR311E	Tek-CARE300III	Single Patient Station, with RJ45 8-pin modular jack	347.89
IR312E	Tek-CARE300III	Dual Patient Station , with two RJ45 8-pin modular jacks	407.97
IR319B	Tek-CARE300III	Single Patient Station, 8-pin DIN	223.66
IR320B	Tek-CARE300III	Dual Patient Station, 8-pin DIN	247.14



## Multi-Purpose Station– OEM


Multi-Purpose Stations provide audible, visual and digital communications, but do not have the controls, LEDs or speakers normally associated with patient, staff or duty stations. Applications may include providing input for code emergency call stations in an ICU or providing input for bath emergency call switches in public bath areas.

- For use with Tek-CARE300III System
- Microprocessor design
- UL® 1069 Listed

Part No.	System	Description	Price
IR300E	Tek-CARE300III	Multi-Purpose Station, 12 VDC	204.75

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.

# Staff Station & Duty Station




### Staff Station – OEM

Staff Pushbutton Station from TekTone. Mounts on a 2-gang box, with a plastic faceplate. Features "call" and "cancel" buttons.

- For use with Tek-CARE300III System
- Call-placed & privacy indicators for ease of use

Part No.	System	Description	Price
IR310C	Tek-CARE300III	Staff Station, pushbutton	245.49



### Duty Station– OEM

Duty station from TekTone. Mounts on a 2-gang box, with a plastic faceplate. Features "call" and "reset" buttons.

- For use with Tek-CARE300III System
- Speaker incorporated
- UL® 1069 Listed

Part No.	System	Description	Price
IR315F	Tek-CARE300III	Duty Station, with speaker	245.49

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.



## Single Patient Station – OEM

The Single Patient Station provides reliable two-way hands-free communication between patient and nurse, and provides for programmable privacy, paging and patient priority from the master station. Call-placed indicator on the panel indicate the call has been placed and the intercom is in use.

- For use with Tek-CARE300III System
- Includes 6' washable TPU plastic pullcord
- UL® 1069 Listed

Part No.	System	Description	Price
IR318B	Tek-CARE300III	Single Patient Station, with call-placed indicator	220.44



JK300A



JK300/2

## Monitor Call Interface Stations – OEM

The Monitor Call Interface from TekTone is for use in locations where a simple interface is needed. Accepts a standard 1/4" plug and registers a "cord out" call. Dummy plug included for use when not connected to a device.

- For use with Tek-CARE300III System
- Wire lead pigtail connections

Part No.	System	Description	Price
JK300A	Tek-CARE300III	Monitor Call Interface Station, 1/4" phone jack	31.73
JK300/2	Tek-CARE300III	Monitor Call Interface Station, 1/4" dual phone jack	50.10

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.





Emergency Pullcord Station – OEM

Emergency Pullcord Station from TekTone features a pull-down style actuator with optional pullcord, and a call placed LED. Mounts on a 1-gang box, with a plastic faceplate.

- For use with Tek-CARE300III System
- Call-placed indicator

Part No.	System	Description	Price
SF340B	Tek-CARE300III	Emergency Pullcord Station	53.44



Code Blue Station– OEM

This TekTone station features a code blue slide activator. A call can be canceled by returning the slide activator to the up position. Includes a call-placed indicator for call assurance. Ideal for ICU/CCU, physical therapy rooms, post-op recovery and more.

- For use with Tek-CARE300III System
- Call-placed indicator
- UL® 1069 Listed

Part No.	System	Description	Price
SF341B	Tek-CARE300III	Code Blue Pullcord Station	63.65

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.



## Staff Presence Station – OEM

The Staff Presence Station provides push button registration of nurse presence (green LED indicator), or aid presence (yellow LED indicator) at remote locations. Can also be used to transfer nurse follower feature from room to room as nursing staff proceeds with rounds. Depressing the nurse presence button automatically transfers indication of all subsequent patient calls to the remote location of the nurse, as well as the master station.

- For use with Tek-CARE300III System
- Plug-in connectors provide for simple wiring
- UL® 1069 Listed

Part No.	System	Description	Price
SF350B	Tek-CARE300III	Staff Presence Station	97.24



## Visual-Only Patient Station – OEM

Visual-Only Patient Stations from TekTone provide routine, bath emergency and code calls to the master station.

- For use with Tek-CARE300III System
- Call-placed indicator
- Low power requirements
- UL® 1069 Listed

Part No.	System	Description	Price
SF380A	Tek-CARE300III	Visual-Only Patient Station	143.50

All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors. V12. Revised 11/28/23.