



## Crest TouchCall™ Call Cords



### Description

The TouchCall provides reliable and versatile nurse call actuation for patients with limited dexterity. The TouchCall features a large, ultra-sensitive touch plate surface that can be easily pressed by a patient hand, foot, elbow, chin, etc.

The tent design, clip and Velcro® attachments allow for multiple positioning possibilities that make it simple and convenient for any patient to use.

### Features

- Available in a variety of plug styles including 1/4", 8-pin DIN RJ45 or other nurse call system styles
- Standard 10 ft length to suit a variety of applications
- Not for use with nurse call systems that require locking devices

### Technical Specifications

<b>Electrical Rating:</b>	28 VAC/DC	<b>Color:</b>	Classic Grey
<b>Electrical Connection:</b>	Various connection types available	<b>Weight:</b>	0.7 lbs
<b>Button Life:</b>	>1,000,000 cycles	<b>Cord Length:</b>	10'
<b>Switch Type:</b>	Normally Open SPST	<b>Dimensions:</b>	2 3/4"(H) x 5 5/8"(W) x 5"(L)
<b>Warranty:</b>	2 Years	<b>Environment:</b>	Operating temperature: 50° - 100°F
<b>Material:</b>	UL 94V-0 Flame Rated, Latex Free PVC		Relative humidity: 30% to 80% (non-condensing)
<b>Standards:</b>	RoHS II Compliant FDA Listed		

**Crest Technical Support: 1-800-328-8908**

Specifications subject to change without notice.  
Visit [CrestHealthcare.com](http://CrestHealthcare.com) for additional technical documents.