



Crest DuraColor™ Momentary Call Cords



Description

Crest's neon green pendant in the DuraCall Call Cord line is the latest option for nurse call visibility and ease of use. Bright and vibrant, the pendant is ideal for those with poor eyesight or for those who have difficulty distinguishing it from bedding.

Features

- High contrast neon color for easy identification
- Prominent color assists patients & staff
- White cable
- FDA Listed

Technical Specifications

Electrical Rating:	24 VAC/DC	Weight:	5.8 oz (7' version)
Electrical Connection:	Various connection types available	Pendant Dimension:	2.75"H x 1.625"W
Button Life:	>1,000,000 cycles	Cord Length:	Various lengths available
Switch Type:	Normally Open SPST	Environment:	Operating temperature: 50° - 100°F
Warranty:	2 Years		Relative humidity: 30% to 80% (non-condensing)
Material:	UL 94V-0 Flame Rated, Latex Free PVC		Oxygen safe: Meets UL 1069 oxygen enriched atmosphere testing
Standards:	FDA Listed RoHS Compliant ETL Recognized to UL 1069		Liquid ingress protection: Meets IEC 60529 IPx4 requirements
Color:	Neon green pendant with red call button		

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.
Visit [CrestHealthcare.com](https://www.CrestHealthcare.com) for additional technical documents.