



Crest EconoCall™ Momentary Call Cords



Description

If you're looking for the most affordable call cord available, the EconoCall delivers. The EconoCall combines a clean, basic design with reliable function at any price any facility can afford. The smooth pendant shape fits comfortably in the hand.

Features

- Our most economically priced call cord
- Repairable button for fast fixes
- Plug styles available to fit most systems
- FDA listed
- RoHS 2 Compliant

Technical Specifications

Electrical Rating:	24 VAC/DC	Weight:	5.8 oz (7' version)
Electrical Connection:	Various connection types available	Pendant Dimension:	2.24"H x 1.10"W
Button Life:	10,000 cycles	Cord Length:	Various lengths available
Switch Type:	Normally Open SPST	Environment:	Operating temperature: 50° - 100°F
Warranty:	1 Year		Relative humidity: 30% to 80% (non-condensing)
Material:	UL 94V-0 Flame Rated, Latex Free PVC		Oxygen safe: Meets UL 1069 oxygen enriched atmosphere testing
Standards:	FDA Listed RoHS Compliant		
Color:	White, grey & beige red call button		

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.
Visit [CrestHealthcare.com](https://www.CrestHealthcare.com) for additional technical documents.