



Crest CleanCall™ Momentary Call Cords



Description

Maintaining a clean environment can be difficult, but the CleanCall helps simplify the job. The completely sealed, membrane style button can be easily wiped down to help prevent the spread of germs, making it a great option for infection control in any facility.

Features

- Smooth membrane is easier to sanitize than regular buttons
- Easy-grasp winged design fits comfortably in the hand
- High-contrast red button with nurse call symbol is highly visible
- FDA listed

Technical Specifications

| | | | |
|-------------------------------|--------------------------------------|---------------------------|---|
| Electrical Rating: | 24 VAC/DC @ 0.1 A | Weight: | 6.08 oz (8' version) |
| Electrical Connection: | Various connection types available | Pendant Dimension: | 2.75"H x 1.625"W |
| Button Life: | >1,000,000 cycles | Cord Length: | Various lengths available |
| Switch Type: | Normally Open SPST | Environment: | Operating temperature: 50° - 100°F |
| Warranty: | 1 Year | | Relative humidity: 30% to 80% (non-condensing) |
| Material: | UL 94V-0 Flame Rated, Latex Free PVC | | Oxygen safe: Meets UL 1069 oxygen enriched atmosphere testing |
| Standards: | FDA Listed | | |
| Color: | White pendant with red call button | | |

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.
Visit [CrestHealthcare.com](https://www.cresthealthcare.com) for additional technical documents.