



Crest Elite Digital™ Pillow Speakers

For XN and EN pillow speaker part numbers

Part number can be located on back of the unit

Features & Benefits

Elite Digital Pillow Speakers are designed to automatically detect any healthcare grade TV using our patented AutoSense™ Technology. They are available with either Direct Channel Access (10-key) or digital channel scroll up/down function.

The pillow speaker has features that provide a superior patient experience, including a lightweight, ergonomic design; crisp, gentle-press buttons; and bright, color-coded functions for easy navigation. They meet UL 1069 requirements, and are safe for use in oxygen enriched environments. The case is UL 94V-0 Flame Rated, Latex Free ABS Plastic.

- FDA Listed
- Superior sound quality that lasts
- Water & contaminant resistant
- Durable case & strain relief
- Available in backlit models
- 2 year warranty

Technical Specifications

| | | | |
|---------------------------|---|---------------------|--|
| Electrical Rating: | 30 VDC at 100 mA | Color: | White |
| Weight: | 18.2 oz | Standards: | ETL Recognized to UL 1069 |
| XN Dimensions: | 3.15"W (tapering to 2.10") x 9"H x 1.37"T | | No. 205 with oxygen enriched environment marking |
| EN Dimensions: | 3.2"W (tapering to 2.10") x 9.35"H x 1.26"T | | RoHS Compliant |
| Cable Length: | 8' | | FDA Listed |
| Material: | UL 94V-0 Flame Rated, Latex Free ABS | | FCC Part 15 subpart B |
| Warranty: | 2 Years | Environment: | ICES-003 |
| Speaker Material: | Moisture-resistant mylar cone, 2 watt | | Operating temperature 50° to 120° F |
| | | | Relative Humidity: 30% to 90% (non-condensing) |

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.

Visit [CrestHealthcare.com](https://cresthealthcare.com) for additional technical documents.