



Crest Elite Infrared™ Pillow Speakers

For XT and ET pillow speaker part numbers

Part number can be located on back of the unit

Features & Benefits

Elite Infrared Pillow Speakers are the only pillow speakers that combine complete TV, DVD and Set Top Box Control. They can be programmed to control any healthcare grade and most consumer grade TVs. They are available with either Direct Channel Access (10-key) or digital channel scroll up/down function.

The pillow speaker has features that provide a superior patient experience, including a lightweight, ergonomic design; crisp, gentle-press buttons; and bright, color-coded functions for easy navigation. They meet UL 1069 requirements, and are safe for use in oxygen enriched environments. The case is UL 94V-0 Flame Rated, Latex Free ABS Plastic.

- FDA Listed
- Superior sound quality that lasts
- Water & contaminant resistant
- Durable case & strain relief
- Available in backlit models
- 2 year warranty

Technical Specifications

Electrical Rating:	30 VDC at 100 mA	Color:	White
Infrared Power:	2 AAA batteries (For Entertainment Functions Only)	Standards:	ETL Recognized to UL 1069 No. 205 with oxygen enriched environment marking RoHS Compliant FDA Listed FCC Part 15 subpart B ICES-003
Weight:	19.6 oz	Environment:	Operating temperature 50° to 120° F Relative Humidity: 30% to 90% (non-condensing)
XT Dimensions:	3.15"W (tapering to 2.10") x 9"H x 1.37"T		
ET Dimensions:	3.2"W (tapering to 2.10") x 9.35"H x 1.26"T		
Cable Length:	8'		
Material:	UL 94V-0 Flame Rated, Latex Free ABS		
Warranty:	2 Years		
Speaker Material:	Moisture-resistant mylar cone, 2 watt		

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.

Visit [CrestHealthcare.com](https://www.cresthealthcare.com) for additional technical documents.