

SHOW YOUR PATRIOTISM AND PRIDE TODAY

Four Star Elite Digital™ & Elite Infrared™ Pillow Speakers

NEW Show your patriotism and pride

Four Star pillow speakers include our high quality features, including:

- Elite Digital with AutoSense™ proprietary technology automatically detects and connects to any hospital grade TV
- Elite Infrared is configurable to control any combination of TVs, set-top boxes and DVD players
- Includes auxiliary nurse call function and two-way intercom support
- Crisp, gentle-press buttons
- Environmental controls
- Easy to repair
- FDA listed
- Infection control

Call today for Crest pillow speaker part number!

LATEX-FREE

OXYGEN SAFE

2
YEAR
WARRANTY

RECOGNIZED
COMPONENT
ETL
Intertek

UL
LISTED

CREST
HEALTHCARE®

Crest Healthcare

Phone: 800-328-8908 Fax: 800-369-9207

CrestHealthcare.com