

Nurse Call System Plug Guide

1/4" Phone Plug
2-conductor plug




1/4" Phone Plug
3-conductor plug



Amplion®
15-pin plug




**Ascom®
Telligence®
Dukane®**
18-pin plug




Critical Alert™
8-pin DIN plug



Dukane®
12-pin plug



Dukane®
12-pin plug



Dukane®
17-pin plug



Ektacom®
8-pin plug



Ektacom®
15-pin plug



Executone®
7-pin plug



**HCE
Sylvania®**
20-pin plug



Hillrom™
20-pin plug



Hillrom™
37-pin plug



Intercall™
7-pin plug



Intercall™
8-pin DIN plug



Jeron®
8-pin DIN plug




Rauland®
8-pin DIN plug




**Rauland®
Responder® 5**
9-pin DIN plug




Rauland®
12-pin plug



RCA®
5-pin plug




RCA®
9-pin plug



**Simplex®
TekTone®
Critical Alert™**
RJ45 8-pin plug




Simplex®
15-pin plug




TekTone®
8-pin DIN plug



Zenith®
6-pin plug



Zettler®
6-pin plug



Zettler®
9-pin plug

