



LATEX-FREE

2
YEAR
WARRANTY

Crest SoftCall™ Momentary Call Cords

Description

SoftCall provides a reliable and affordable option when only nurse call and light control is needed. The waterproof pendant fits comfortably in the hand and can be easily wiped down for infection control. The high-contrast color makes the SoftCall pendant easy-to-find in bedding.

Features

- High contrast blue color for easy identification
- Soft silicone pendant
- 8' white cable
- FDA Listed
- IP47 sealed & waterproof

Technical Specifications

Electrical Rating:	24 VAC/DC	Weight:	9.6 oz
Electrical Connection:	Various connection types available	Pendant Dimension:	4.3"L x 1.0"H x 1.8"W
Button Life:	>50,000+ cycles	Cord Length:	8'
Switch Type:	Normally Open SPST	Environment:	Operating temperature: 50° - 100°F
Warranty:	2 Years		Relative humidity: 30% to 80% (non-condensing)
Material:	Silicone, Latex Free		Liquid ingress protection: Meets IEC 60529 IP47 requirements
Standards:	FDA Listed RoHS Compliant		
Color:	Blue pendant with red call button and yellow light buttons		

Crest Technical Support: **1-800-328-8908**

Specifications subject to change without notice.
Visit [CrestHealthcare.com](https://www.cresthealthcare.com) for additional technical documents.