

Features

- Enables dual mounting of SF-series stations on the Tek-CARE120 and 160 systems.

Description

The IH122K mount kit allows two SF-series Tek-CARE120 or 160 stations to be mounted in a single housing. Any combination of two SF121, SF122, and SF123 stations may be installed in the IH122K to satisfy the needs of the facility.

The plastic bezel snaps onto the included mounting plate, leaving no visible screw holes.

Listings



The Tek-CARE120 and 160 nurse call systems are listed by Underwriters Laboratories, Inc. to the UL®1069 seventh edition standard for Hospital Signaling and Nurse Call Equipment, the UL®2560 standard for Emergency Call Systems, and cUL® Listed to CSA C22.2 No. 205. UL® File Number E71599.



RoHS Compliant.



IH122K

Specifications

- Dimensions: 4.58"H × 6.74"W × 1.03"D (116 mm × 171 mm × 26 mm)
- Environment: 0°C–60°C
- Construction: Flame-retardant (UL® 94V-0) plastic faceplate.
- Mounting: Attach included back plate to a standard single- or double-gang ring, or single- or double-gang ring and box.

Required Components

- SF121 Single Patient Station
- SF122 Dual Patient Station
- SF123 Two-Button Pull-Cord Station
- Included with the SF125 Duty Station