

Premium Raised Toilet Seat with Removable Arms

HCPCS: E0244



FIGURE A



FIGURE B

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Allows individual to use existing toilet seat and lid
- Tool-free arms are 19.75" apart and can be removed or added as needed
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Lightweight and portable
- Easily attaches to toilet bowl
- No tools required for installation
- Easy to clean
- Two models accommodate standard and elongated toilets
- Made using recycled materials

PRODUCTS

Description

12402	Rizer,RegularToiletsKD,w/RemArm,1/cs,RTL
12403	Rizer,ElongatedToilet,w/RemArms,1/cs,RTL

SPECIFICATIONS

Description

Carton Shipping Weight	4.5 lbs.
Seat Dimensions	13.75" (W) x 17" (D) x 3.5" (H) 13.75" (W) x 19.25" (D) x 3.5" (H)
Warranty	Limited Lifetime
Weight	3.5 lbs.
Weight Capacity	300 lbs
Width Between Arms	19.75"