


# Prevent Falls & Wandering

## CordLess Pads & Mats For Bed & Chair Exit Alarms



Works with:  
TL-2100G



Works with:  
433-EC  
433-CMU-40/60

### No Cords Between Pads & Monitors

For a quieter environment, use a CordLess Pad with a CordLess Monitor. This allows you to remove the alarm noise from the room. The Pad is placed under resident, either on a chair or bed, and the floor mat can be placed either by the bed or by the door. When a resident gets up from a pad, or steps on a floor mat, a wireless signal is sent to the monitor to notify caregivers.

Works with:



CordLess  
Monitor  
TL-2100G

Works with:



Central  
Monitor Unit  
433-CMU-40/60



CordLess  
Alarm  
433-EC

MFR#	SMART PRODUCT DESCRIPTION
GBT-WI	CordLess Bed Pad 20"x30" (1 year)
GBT-RI	CordLess Bed Pad 10"x30" (1 year)
GBT-90, GBT-45	CordLess Bed Pad 10"x30" (90 day, 45 day)
GCT-WI	CordLess Chair Pad 10"x15" (1 year)
GCT-90, GCT-45	CordLess Chair Pad 10"x15" (90 day, 45 day)
FMT-05C	CordLess Floor Mat 24"x36" Gray (1 year)
FMT-07C	CordLess Floor Mat 24"x48" Gray (1 year)
GWC-02	CordLess Fall Mattress 70.5"x29.5" x1.5" (90 day)
LM-01C	CordLess Landing Mat 24"x71"x1" Gray (1 year)

MFR#	SMART PRODUCT DESCRIPTION
PTB-WI	Economy CordLess Bed Pad 20"x30" (1 year)
PTB-RI	Economy CordLess Bed Pad 10"x30" (1 year)
PTB-90R, PTB-45R	Economy CordLess Bed Pad 10"x30" (90 day, 45 day)
PTC-WI	Economy CordLess Chair Pad 10"x15" (1 year)
PTC-90, PTC-45	Economy CordLess Chair Pad 10"x15" (90 day, 45 day)
PTFM-05C	CordLess Floor Mat 24"x36" Gray (1 year)
PTFM-07C	CordLess Long Floor Mat 24"x48" Gray (1 year)
PTWC-02	CordLess Pressure Fall Mattress 70.5"x29.5"x1.5" (1 year)
PTLM-01C	CordLess Weight Landing Mat 24"x71"x1" (1 year)

(800) 650-3637



SmartCaregiver.com