

Gel "E" 2" General Use Gel/Foam Wheelchair Cushion

HCPCS: E2601




FIGURE A



FIGURE B

- The Medicare approved General Use 2" gel-foam wheelchair cushion will provide comfort and support for the wheelchair bound patient
- High-Density foam encapsulates a viscous gel bladder that provides proper immersion and envelopment for pressure redistribution
- The cover is fluid resistant low-shear nylon zippered and removable with a durable vinyl bottom
- Carry handle included for convenient transport
- Product Dimensions: 18" (W) X 16" (D) X 2" (H)
- Meets flammability standard TB117

PRODUCTS

Description

| | |
|-------|---|
| 14888 | W/C Cushion,2in Gel/Foam,18x16x2in,1/cs |
|-------|---|

SPECIFICATIONS

Description

| | |
|------------------------------------|---|
| Cover Material | Fluid-Resistant Nylon Waterproof Vinyl Bottom |
| Warranty | 1 Year Limited |
| Weight Capacity | 250 lbs |
| Weight Capacity (13046 13046-HD) | 250 lbs. |