



**The Aurora from SMS Technologies is the most powerful handheld, portable UVC disinfection device on the market.**

The Aurora delivers  $>10,000\mu\text{W}/\text{cm}^2$  UV exposure at 8"

The powerful design deactivates most common viruses and bacteria in around one second of exposure.

**Kill times for common pathogens**

4 log kill (99.99% deactivation),  
Aurora at 8" distance


	Seconds
Coronavirus	<1
Influenza	1
MRSA	1
Hepatitis	2
Norovirus	3
CDiff	8

**"This thing is a cannon!"**

**"In lab testing, all the plates that were set up and exposed to the UV light failed to grow. It works as advertised!"**

*– Hospital lab manager*

**CREST**  
HEALTHCARE®

 800.328.8908

 800.369.9207

 CrestHealthcare.com

**SMS**  
TECHNOLOGIES