



Crest Break-A-Way Communication Cables

Description

Crest Break-A-Way Bed Communication Cables utilize a pinless 36-pin Centronix connector two feet from the bed plug end. If the bed is moved without disconnecting, the cable automatically breaks away at the connector, saving the cable and wall station from damage.

These sturdy cables are available for Hill-Rom® or Stryker® beds with in-line or right angle plugs.

Features

- Saves stations from damage
- Economical & durable
- Break-A-Way located 2' from bed end
- 8' beige cord

Technical Specifications

Electrical Rating:	24 VAC/DC	Color:	Beige
Electrical Connection:	37-pin to 36-pin Centronix	Weight:	1.4 lbs
Warranty:	1 Year	Cord Length:	8'
Material:	UL 94V-0 Flame Rated, PVC	Environment:	Operating temperature: 50° - 100°F Relative humidity: 30% to 80% (non-condensing)
Standards:	RoHS II Compliant FDA Listed		

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.
Visit CrestHealthcare.com for additional technical documents.