

GripProtect[®]

Procedure Earloop Mask



Comfortable 3-level Protection For Your Face

GripProtect Procedure Earloop Masks are designed using 3 protective layers for filtration against bacteria, particulates and fluids.

These masks meet ASTM level 1 testing standards for light protection.

Features:

- Elastic earloops & an aluminum nose strip permit a close, comfortable fit
- Latex-free
- Single-use only, disposable
- Unisex, fits most faces
- Blue color, 50 masks per box



For more information visit:
www.baymedicalcompany.com
1-800-977-7888

GripProtect®

Procedure Earloop Mask

- Meets ASTM Level 1 standards
- Safe, comfortable and cost effective
- Comfortable elastic earloops
- Disposable, single use
- Aluminum strip for a close fit over the nose
- Easy to breathe when wearing mask
- Latex-free
- 3-Ply layered design
- 3 Year shelf life



Specifications	
Material—Outer layer	25 GSM (grams per square meter) Spun-bonded nonwoven materials
Material—Middle layer	25 GSM Melt blown; Melt blown nonwoven materials (BFE: 99%)
Material—Inner layer	25 GSM Spunbonded nonwoven materials
Performance standards	ASTM Level 1
Size of each mask (flat)	17.5 cm wide x 9.5 cm height
Quantity / Box	50
Quantity / Case	20 Boxes/Case; 1,000 Masks/Case
Case size	20.5" x 15" x 7.1" (52 x 38 x 18 cm)
Box size	7.1" x 3.9" x 3.1" (18 x 10 x 8 cm)
Part Numbers	Case: MASK3550 Box: MASK3550-BX
Country of origin	China
UPC for case	10850021591127
UPC for box	850021591120



For more information visit:
www.baymedicalcompany.com
 1-800-977-7888