

Features

- Enables single mounting of SF-series stations on the Tek-CARE Nurse Call Systems.

Description

The IH121K mounting kit allows a single SF-series Tek-CARE120 or 160 station to be mounted in a single housing. An IR160, SF121, SF122, or SF123 station may be installed in the IH121K to satisfy the needs of the facility.

The plastic bezel snaps onto the included mounting plate, leaving no visible screw holes. The pull-switch slider and pull-string are included in the kit. The IH121K is packaged with a locking screw for added security.



IH121K Single Station Mounting Kit
(shown without pull-string and slider)

Specifications

Dimensions: 4.58"H × 3.25"W × 1.03"D
(116 mm × 83 mm × 26 mm)

Environment: 0°C–60°C

Construction: Flame-retardant (UL® 94V-0) plastic.

Mounting: Attach included back plate to a standard single-gang ring, or single-gang ring and box.

Required Components

IR160 Audio Station
SF121 Single Patient Station
SF122 Dual Patient Station
SF123 Two-Button Pull-Cord Station