

LifeGlider™ Mobility Device



116141

LifeGlider™ Mobility Device

The LifeGlider is a new mobility device that allows upright, hands-free mobility. It is designed to reduce the fear of falling and fatigue sometimes associated with struggling to walk. This patented assistive device helps users maintain balance, prevents falls and move independently. It can support all, some or none of the user's weight.

- Useful in therapy settings & for ongoing personal use
- Lightweight & portable
- Can be used indoors & outdoors in any ADA compliant environment

200
LB WEIGHT
CAPACITY

275
LB WEIGHT
CAPACITY

1
YEAR WARRANTY ON
PARTS & LABOR

1
YEAR WARRANTY
ON FRAME

Part No.	Description	Price
* 116261	LifeGlider, small, black, 200 lb	695.95
* 116141	LifeGlider, standard, black, 275 lb	695.95
* 116142	LifeGlider, tall, black, 275 lb	695.95

SPECIFICATIONS

	116261	116141	116142
Front Width:	26"	26"	26"
Caster Size:	8"	8"	8"
Seat-to-Floor Height:	27" – 32"	30" – 40"	33" – 40"
User's Height:	4'6" – 5'2"	5'3" – 5'9"	5'10" – 6'4"
Minimum Inseam:	24" – 27"	28" – 32"	33" – 38"
Weight:	17 lb	19 lb	21 lb

CREST
HEALTHCARE®

Crest Healthcare

Phone: 800-328-8908 Fax: 800-369-9207
CrestHealthcare.com

©2020 Crest Healthcare. *Shipped direct from the supplier to you.
All listed prices and products are subject to change without notice.
Crest is not responsible for photographic or typographical errors.