

Crest Remote Call-IR™

XP441EOB600XXX2

XP441XOB600XXX2

Remote Call-IR™ is a new nurse call device that gives residents and staff the convenience of a TV remote control combined with a call cord. Choose either 10-key direct channel access or large TV buttons. The large button model allows the resident to easily control TV power, channel and volume with an IR signal. It's tethered to the wall, by the nurse call cable always staying within easy reach. Remote Call-IR eliminates the extra buttons on remotes that can frustrate residents. Can be programmed for use with multiple TV brands like a universal remote.

For safety, the nurse call function requires no battery power and will function reliably even if battery life runs out. The TV controls are powered by 2 AA batteries (not included) that will last about as long as typical remote control batteries. Remote Call-IR includes a 10' cable and a 2-conductor 1/4" phone plug to connect to the nurse call systems. Remote Call-IR was designed to comply with UL 1069 requirements and is safe for use in oxygen-enriched environments.

Specifications

ELECTRICAL**Current and****Voltage Rating:**

1.5A @ 24V

OXYGEN SAFETY:

Meets UL 1069 Requirements

LATEX INFORMATION:

This product is latex free

MECHANICAL**Button Life:**

>1,000,000 cycles

Durability:

100 Drops on Concrete from 4' height

SIZE MEASUREMENTS**Height:**

8.9"

Width:

3.5"

Thickness:

1.82"

WARRANTY:

2 Year

**10-Key Model****Big Button Model**