



Expand Your Wireless Care

Crest Wireless System Extender, for multiple transmitters

The Wireless System Extender, for multiple transmitters, allows the option of expanding existing nurse call stations without re-wiring. The wireless receiver connects to a universal 1/4" plug station and can be used with a wearable pendant, wireless Fall Fighter monitor, wireless smoke detector, or a wireless nurse call station. Program up to 4 wireless nurse call devices to the receiver. Pendant battery life is 1-2 years, depending on usage. All wireless devices include batteries. Replacement batteries are also sold separately.

- **No need to damage walls to add stations to resident room**
- **Program up to 4 wireless nurse call products to 1 receiver**
- **Easy to install & use**
- **Repairable by Crest Repair Technicians**

Part No.	Description	Price
116121 [†]	Wireless System Extender, for multiple transmitters, with power supply and 8' 1/4" plug	298.73
*EN1223S	Wireless Pendant, water resistant, with neck cord and belt clip	93.79
116122	Wireless Fall Fighter Monitor	167.43
116123	Wireless Station, momentary pull cord	188.95
116124	Wireless Station, momentary 1/4" jack	162.02
EN1244	Wireless Smoke Detector	192.99
BAT609	Replacement Battery, wireless pendant	4.14
114820	Replacement Battery, wireless stations	4.95



Additional wireless system details available
View entire selection details online at CrestHealthcare.com

CREST
HEALTHCARE®

 800.328.8908

 800.369.9207

 CrestHealthcare.com