



Crest DuraCable™ Bed Communication Cable

Description

Crest DuraCables are the affordable and durable option for beds that have communication functions built into the side rail or a bed exit alarm feature.

Hill-Rom® cables feature one male plug end and one female plug end which connect to 37-pin wall outlets.

Stryker® cables feature 37-pin male plugs on both ends.

Features

- In-line or right angle plug options
- Molded plug ends with strain relief
- Available for Hill-Rom® & Stryker® beds

Technical Specifications

Electrical Rating:	24 VAC/DC	Color:	Beige
Electrical Connection:	Various connection types available	Weight:	1.2 lbs
Warranty:	1 Year	Cord Length:	8'
Material:	UL 94V-0 Flame Rated, PVC	Environment:	Operating temperature: 50° - 100°F Relative humidity: 30% to 80% (non-condensing)
Standards:	RoHS II Compliant FDA Listed		

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.
Visit CrestHealthcare.com for additional technical documents.