



MJM International - The Leader In Designing And Manufacturing Of PVC Medical Products


www.mjminternational.com

New And Innovated Gurney Duo Sling Design

Advantages

- PADS ARE SEWN ONTO THE MESH SLING AND STAY IN PLACE WHILE GURNEY IS IN USE
- SAFER OPTION FOR END USER AND AIDS
- IMPROVED DRAINAGE, WATER WILL NO LONGER PUDDLE
- INCREASED PATIENT COMFORT
- COST EFFECTIVE, SHIPS IN A SMALLER BOX THUS SAVING ON FREIGHT
- EASY TO CLEAN AND MAINTAIN
- DESIGNED TO LAST LONGER THAN STANDARD FOAM PAD

Mesh Fabric
Between Padding
For Maximum
Drainage



Closed Cell Foam
Pads Are Sewn
Onto Gurney Sling
For Superior Comfort



For More Information Call MJM International Today!

800-270-6199 Fax: 956-928-9535 www.mjminternational.com sales@mjminternational.com