

2500KL
Digital Wheelchair Ramp



2500CKL
Digital Wheelchair Ramp Scale
with Fold Away Seat



Advantages

- ✓ **EMR Connectivity**
Scale interfaces with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Easy Portability**
The scale has two wheels that provide easy mobility and allow for convenient placement
- ✓ **Swivel Display Head**
For the convenience of the administrator, the display head rotates and pivots allowing line-of-sight viewing and information privacy
- ✓ **On-Screen Help Menu**
Guides users through operating functions and provides easy access to multiple settings
- ✓ **Everlock® (patent pending)**
Option that allows users to permanently lock the scale in either LB or KG to improve safety
- ✓ **Interactive Display**
High-contrast color screen provides a user-friendly interface
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost



2500CKL / 2500CKG

2500KL

2500CKL

Specifications



- Capacity: 1,000 lb / 454 kg
- Seat Capacity: 419 lb / 190 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
29 3/8" x 28 1/4" x 2 3/8" / 746 mm x 718 mm x 60 mm
- Seat Size for 2500CKL, 2500CKG or CHAIRACC (w x d x h):
18" x 15" x 18 1/2" / 457 mm x 381 mm x 470 mm
- Ramp Size (w x d x h): 28 1/4" x 9 7/8" / 714 mm x 248 mm
- Product Footprint (w x d x h):
40 1/4" x 34 1/2" x 45 1/2" / 1022 mm x 876 mm x 1156 mm
- Product Weight: 121 lb / 55 kg
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Wheels: 2
- Power Source: 100-240V adapter included, 6 D batteries (not included)
- Connectivity: USB, optional Pelstar wireless technology
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off, Time, Audible/Mute Sound Option, Everlock®
- 2 year limited warranty

Model	Description
2500KL	Digital wheelchair ramp scale, pounds and kilograms
2500KL-BT	2500KL with Pelstar wireless technology
2500KG	Digital wheelchair ramp scale, kilograms only
2500KG-BT	2500KG with Pelstar wireless technology
2500CKL	Digital wheelchair ramp scale with fold away seat, lb & kg
2500CKG	Digital wheelchair ramp scale with fold away seat, kilograms only
SS-2500KL	ScaleSurance 2 year extended warranty for 2500KL, 2500KL-BT, 2500KG, & 2500KG-BT*
SS-2500CKL	ScaleSurance 2 year extended warranty for 2500CKL & 2500CKG
CHAIRACC	Optional fold away chair attachment for 2500KL, 2500KL-BT, 2500KG, & 2500KG-BT
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
B2500RAMP	Optional secondary ramp for 2500 models
PELSTARONE-P	Wireless thermal printer