

# Crest Pillow Speaker Tester



**1**  
YEAR  
WARRANTY

Call Crest's Technical Support Team To  
Discuss Your Facility Needs

Test functionality of pillow speakers as well as bed controls and call cords with this compact and portable tester.

- Tester works with analog, digital & infrared pillow speakers
- Compact & portable
- Operates with 12 VDC or 9 V battery (both included)
- Modules to fit nurse call systems sold separately

Part No.	Description	Price
115642	Pillow Speaker Tester, white, not compatible with Rauland® Responder® 5	238.95
115643	Pillow Speaker Tester, white, compatible with Rauland® Responder® 5	260.95
CALL	Standard Module, SST	130.95
CALL	Premium Module, SST	151.95

## SPECIFICATIONS

<b>Length:</b>	5 ¾"
<b>Width:</b>	4 ⅝"
<b>Height:</b>	1 ¾"
<b>Power Requirements:</b>	9 V Battery or 12 VDC adapter (both included)

**CREST**  
HEALTHCARE®

**Crest Healthcare**

Phone: 800-328-8908 Fax: 800-369-9207

[CrestHealthcare.com](http://CrestHealthcare.com)

©2019 Crest Healthcare. All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors.