

# Broken Bed Cables? Call CREST!

## BLACK JACK

### A Magnetic Connection

BlackJack Magnetic Breakaway Bed Cables are designed to be universally applicable to most bed types and nurse call systems. The magnetic "brick" adapter is simply screwed into existing wall or bed jacks or onto the end of a standard adapter cable matching the pin-out on the wall. The magnetic breakaway side of the cable attaches to the magnetic brick adapter and guarantees 100% connectivity. Providing a quick disconnect from the wall or bed side, cables automatically break away from the jack to prevent damage if bed is moved without disconnecting. Call Crest to learn more about this magnetic cable solution for your facility.

- Breakaway location can be placed within 2' of the bed or wall
- Works with Hillrom™, Linet®, Stryker® (requires 10043 gender changer), CHG® & 37-pin port beds
- Gender changer & dummy plug sold separately online



MAGNETIC BREAKAWAY "PIGTAIL" SOLUTION



MAGNETIC DUMMY PLUG  
116110

### FEWER BROKEN CABLES = MORE SAVINGS!!

**New Cable Cost:** \$35

**New Station Cost:** \$400

**Lost Revenue from out-of-service room:** \$1500

**TOTAL SAVINGS: \$1935<sup>+</sup>**


<sup>+</sup>Total savings per bed based on estimated cost of new cable, station, and 1-day of lost room reimbursement revenue. Actual savings to facility may vary by cable cost, station cost and daily room charge.




Contact Crest Customer Service at **1-800-328-8908** to get started with your order today!

Part No.	Description	List Price
<a href="#">*116673</a>	BlackJack 37-pin Right Angle Magnetic Cable with Magnet	<b>169.17</b>
<a href="#">*116672</a>	BlackJack 37-pin Magnetic Cable	<b>169.17</b>
<a href="#">*116110</a>	BlackJack Magnetic Dummy Plug	<b>61.79</b>



 800.328.8908

 800.369.9207

 CrestHealthcare.com

©2023 Crest Healthcare. \*Shipped direct from the supplier to you. All listed prices and products are subject to change without notice. Crest is not responsible for photographic or typographical errors.