



Wireless Fall Fighter® Monitor



RoHS
Compliant

LATEX-FREE

1
YEAR
WARRANTY

Description

The Wireless Fall Fighter monitor has the robust features you need for chair or bed applications at an affordable price. The alarm automatically resets when a patient returns to the pad or is quickly silenced with the reset button. A battery cover with screws prevents tampering. The front LED flashes blue to indicate low battery and red to indicate active status and alarm. The monitor includes a nurse call interface jack, AC adapter input and volume options to silence alarm in room when connected to a nurse call system. Uses two 3 V batteries for mobile use, or optional plug in adapter for stationary use. Monitor protector with fastening strap sold separately.

Features

- Wireless transmitter ready
- “Always on” monitoring
- Chime, warble & buzz alarm tones
- 0 – 8 second delay reduces false alarms

Technical Specifications

Electrical Rating:	Two 3 V (CR123 Type Lithium) batteries or DC 6 V adapter (not included)	Color:	White
Volume Range:	85 – 105 dB, mutable	Weight:	5.3 oz
Material:	ABS	Dimensions:	125 x 76 x 35 mm
Standards:	RoHS Compliant FDA Listed FCC Part 15 subpart B UDI Labeled	Warranty:	1 year
		Status LED Functionality:	Blue - low battery Red - active status and alarm

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.
Visit CrestHealthcare.com for additional technical documents.