

Wedge Cushion with Stretch Cover

HCPCS:



FIGURE A

- Wedge design helps prevent patient from sliding forward
- Firm foundation reduces pressure caused by sling-seat hammocks
- The top cover is fluid resistant, low shear stretch fabric
- The bottom cover is waterproof
- High-density foam increases stability and comfort

PRODUCTS

Description

| | |
|------|--|
| 8134 | WedgeWheelchairCushion18x16x3-1 1ea,4/CS |
|------|--|

SPECIFICATIONS

Description

| | |
|-----------------|---|
| Cover Material | Fluid-Resistant Stretch Top Waterproof Vinyl Bottom |
| Warranty | 1 Year Limited |
| Weight Capacity | 250 lbs |