

Gel "E" Skin Protection 3" Gel/Foam Wheelchair Cushion

HCPCS: E2603




FIGURE A



FIGURE B

- Pressure redistributing cushion designed to assist in the prevention, treatment and management of pressure injuries
- Viscous gel bladder is encased in polyurethane foam shell and provides optimal pressure redistribution, support and comfort
- Cover top is a urethane-coated nylon that is low-shear, vapor permeable and water-resistant
- Cover base is waterproof vinyl for durability
- Carry handle included for convenient transport
- Meets flammability standard TB117

PRODUCTS	Description
14878	Gel "E" W/C Cushion, 18x18x3 ,1/ea, 4/cv
14886	\\I "E" W/C Cushion, 18x16x3, 1/cv, 4/ct
14891	Gel "E" W/C Cushion, 16x16x3, 1/bx, 4/cv
14894	Gel "E" W/C Cushion, 20x16x3, 1/cv, 4/ct
14902	Gel "E" W/C Cushion, 22x18x3, 1/cv, 4/ct

SPECIFICATIONS	Description
Cover Material	Fluid Resistant Nylon Top Waterproof Vinyl Bottom
Warranty	18 Months Limited
Weight Capacity	500 lbs
Weight Capacity (13046 13046-HD)	Up to 20" 275 lbs. 22" and up 500 lbs.