

Gel-U-Seat Skin Protection 3" Gel/Foam Wheelchair Cushion

HCPCS: E2603



FIGURE A

- The Gel-U-Seat skin protection cushion has a high density foam shell surrounding a uniquely designed gel bladder
- Dual chambered gel bladder provides proper immersion and envelopment for pressure redistribution
- Cover is water resistant and low shear stretch to reduce heat and moisture build-up
- Convenient 12" tie straps secure the cushion to the wheelchair
- Meets flammability standard TB117

PRODUCTS

	Description
8046-1	Gel-U-Seat W/C Cush, 16x16x3, 1ea,4/CS
8046-2	Gel-U-Seat W/C Cush, 18x16x3, 1ea, 4/CS
8046-3	Gel-U-Seat W/C Cush, 20x16x3, 1ea, 4/CS
8046-4	Gel-U-Seat W/C Cush, 20x18x3, 1ea, 4/CS
8046-5	Gel-U-Seat W/C Cush, 22x18x3, 1ea,4/CS
8046-7	Gel-U-Seat W/C Cush, 24x18x3, 1ea, 4/CS

SPECIFICATIONS

	Description
Cover Material	Fluid Resistant Stretch Top Waterproof Vinyl Bottom
Warranty	1 Year Limited
Weight Capacity	Up to 20" 275 lbs, 22" and up 500 lbs