

# Stick-On Bags

## “Stick It On...Hang It Anywhere”

- Red Infectious Waste Bags and Beige Regular Waste Disposal Bags
- Red Infectious Waste Bags feature the printed universal biohazard symbol and bilingual biohazard text
- Stick-On Bags have two adhesive strips:
  - one for sticking bag on counter, wall, etc.
  - one to seal the bag when it is full
- Adhesive strips adhere to most clean, dry surfaces
- 2 mil bag thickness
- 2 sizes to accommodate all your waste disposal needs



	Item Number	Description	Quantity
A	MRWB142316	9" x 10" Red Biohazard Waste Stick-On Bag, 1.4 qt, 2 mil	100 ea/bx
	MRWB142324	12" x 14" Red Biohazard Waste Stick-On Bag, 2.6 qt, 2 mil	100 ea/bx
B	MBWB142325	9" x 10" Beige Standard Waste Stick-On Bag, 1.4 qt, 2 mil	100 ea/bx

Product appearance and specifications are current at the time of printing, subject to change without notice. Availability of certain products may be limited by federal, state, provincial or local licensing requirements.



21875 Grenada Avenue, Lakeville, MN 55044  
 toll free: 800-347-9023 • ph: 952-469-9400 • fax: 952-469-9490  
 customercare@unimedcorp.com • www.unimedcorp.com