

NITRO™

3-WHEEL ROLLATOR MOBILITY



Ergonomic handgrips manufactured with Drive PreserveTech™, an antimicrobial* agent that helps to minimize the growth of bacteria, mold and mildew on treated surfaces



Push down on the handbrake to engage parking brake feature



Adjust the handle height with the push of a button



Lightweight aluminum frame with built-in brake cable for extra safety and protection



Large, removable storage pouch is perfect for transporting items to and from



22" width is ideal for tight to navigate spaces such as doorways, hallways and compact living quarters



Dual front fork design enhances turning radius while providing added safety



Large 9" front wheel allows for optimal steering and maneuverability over most types of terrain



Park and store the rollator anywhere with easy side-to-side fold

ITEM #	PRODUCT DESCRIPTION	UPC	UOM
RTL10266TWHL	Nitro 3-Wheel Rollator, Red	822383006314	1/ea

The Nitro 3-Wheel Rollator is the newest addition to our premium line of Nitro rollators. Featuring a hidden inner brake cable, the functional and modern Nitro 3 is designed with safety in mind while indoors or outdoors. The Nitro 3 easily navigates narrow spaces and folds flat while in the standing position so it can roll, park and store anywhere. Also included is a large, removable pouch to keep personal items close. The Nitro 3 comes pre-assembled, so you can unbox, unfold and get your customers going!

*Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

SPECIFICATIONS

Frame	Aluminum
Width Between Arms	12"
Width / Depth	22" (W) x 29" (D)
Handle Height Min/Max	32.50" (min) 37.50" (max)
Wheel size	9"
Product Weight	12.70 lbs.
Color	Red
Carton Dimensions	32" (L) x 11" (W) x 34" (H)
Carton Shipping Weight	17 lbs.
Weight Capacity	300 lbs.
Warranty	Limited Lifetime, Brake and Handle Housing 5 Years

For more information or to place an order, please contact customer service today.

877-224-0946  www.drivemedical.com

©2020 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies.