



VIPER™ PLUS GT WITH UNIVERSAL ARMRESTS

THE BEST SEAT IN THE HOUSE

VERSATILE WHEELCHAIR

NEW

HEAVY-GAUGE, FLAME-RETARDANT NYLON UPHOLSTERY WITH MODESTY FLAP



HEIGHT ADJUSTMENT 8.5-11.5"



ADJUST DESK ARM TO FULL ARM



ALL-ALUMINUM LIGHTWEIGHT FRAME

OVER-CENTER CROSS BRACE ELIMINATES NEED FOR SEAT GUIDES

TOOL-FREE, ADJUSTABLE ELEVATING LEG RESTS WITH EASY-TO-ACCESS RELEASE HANDLES

URETHANE REAR TIRES & CASTERS

PRECISION SEALED FRONT & REAR WHEEL BEARINGS

PUSH-TO-LOCK WHEEL LOCKS

ADJUSTABLE ANGLE & 3 HEIGHT CASTER FORKS



2" Deeper



16 | 18 | 20 | 22"

Four seat widths available



Dual axle allows for multiple seat-to-floor height positions



5-20°

17-19"

Adjustable back angle & height



Adjustable anti-tipper with flip-up wheels



24"

Quick-release, composite, lightweight, maintenance-free Mag-style wheels

ITEM #	PRODUCT DESCRIPTION	UPC	UOM
PLA416FBUARAD-SF	16" Flip Back Universal Arm, Swing-away Footrests	822383005256	1/ea
PLA416FBUARAD-ELR	16" Flip Back Universal Arm, Elevating Legrests	822383005263	1/ea
PLA418FBUARAD-SF	18" Flip Back Universal Arm, Swing-away Footrests	822383005270	1/ea
PLA418FBUARAD-ELR	18" Flip Back Universal Arm, Elevating Legrests	822383005287	1/ea
PLA420FBUARAD-SF	20" Flip Back Universal Arm, Swing-away Footrests	822383005294	1/ea
PLA420FBUARAD-ELR	20" Flip Back Universal Arm, Elevating Legrests	822383005300	1/ea
PLA422FBUARAD-SF	22" Flip Back Universal Arm, Swing-away Footrests	822383005317	1/ea
PLA422FBUARAD-ELR	22" Flip Back Universal Arm, Elevating Legrests	822383005324	1/ea

PRODUCT SPECS	16"	18"	20"	22"
Width Between Arms	16"	18"	20"	22"
Overall Width	24"	26"	28"	30"
Seat Dimension	16" x 16"-18"	18" x 16"-18"	20" x 16"-18"	22" x 18"
Seat-to-Floor Height	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"	17.5"-19.5"
Product Weight	34.5 lb.	35.25 lb.	36 lb.	36 lb.
Color	Deep Green	Deep Green	Deep Green	Deep Green
Weight Capacity	300 lb.	300 lb.	300 lb.	350 lb.
Warranty	Limited Lifetime*	Limited Lifetime*	Limited Lifetime*	Limited Lifetime*

BOX SPECS	16"	18"	20"	22"
Carton Box Dimensions (WxLxH)	11.4" x 31.11" 37"	11.4" x 31.11" 37"	11.4" x 31.11" 37"	11.4" x 31.11" 37"
Carton Box Weight (lb.)	43 48.3	43.3 48.6	43.6 48.9	44.5 49.8

* Refer to instruction guide or website for warranty details.