

# GripProtect® Vinyl Exam Gloves

*Frequent glove changes call for the value and comfort of GripProtect*



GripProtect Vinyl Exam Gloves are high quality, disposable gloves for when simple hand protection and reducing transmission of infectious agents are of primary importance. These gloves are made from durable, synthetic polyvinyl chloride (PVC) and offer a great option in settings where frequent changing of gloves is required.

Classified as Medical Exam Grade, the gloves are manufactured to meet the ASTM testing standards to protect healthcare professionals and their patients from infection. They are recommended for health, lab, food preparation and other uses that do not involve exposure to strong chemicals.

- Medical exam grade gloves
- Powder-free and latex-free
- Clear color, ambidextrous
- Easy to don with smooth inner surface



For more information visit:  
[www.baymedicalcompany.com](http://www.baymedicalcompany.com)  
1-800-977-7888

# GripProtect® Vinyl Exam Gloves

- Medical exam grade
- Latex-Free, Powder-Free
- Non-sterile
- Ambidextrous
- Clear Color
- Easy to don with smooth inner surface



Part Number	Size	# Gloves	Weight / Glove (g)	Width at Palm (cm)
GP1401	Small	100/Box, 1000/Case	4.4 ± 0.3	8.5 ± 3
GP1402	Medium	100/Box, 1000/Case	4.9 ± 0.3	9.5 ± 3
GP1403	Large	100/Box, 1000/Case	5.3 ± 0.3	10.5 ± 3
GP1404	Extra-Large	100/Box, 1000/Case	5.5 ± 0.3	11.5 ± 3

## Specifications, Physical Properties, & Quality Standards

Characteristic	Standard	Test Method
Tensile Strength	≥ 12 MPa	ASTM D5250
Elongation before aging	≥ 340%	ASTM D5250
Thickness at middle finger	0.06 mm (minimum)	ASTM D5250
Thickness at palm	0.07 mm (minimum)	ASTM D5250
Barrier properties	AQL .1.5	ASTM D5250
Length	245 mm	ASTM D5250



For more information visit:  
[www.baymedicalcompany.com](http://www.baymedicalcompany.com)  
 1-800-977-7888