

## UVC Handcart

# Maximize Power. Maximize Mobility.

Adopting UVC Cleaning Systems disinfection technology will complement your existing workflow processes without requiring additional full-time employees. Our products allow your current staff to work smarter, not harder with easy to operate, tested and proven, medical grade, portable UVC Disinfection devices.

## UV EFFECT of 254nm WAVELENGTH

Photochemical destruction of the microorganisms' DNA by absorption of UV-C light which shuts down the organism's reproductive processes rendering it non-pathogenic.

### APPLICATIONS

- Hospitals, Clinics, Medical
- Schools, Commercial
- Home, Residential
- Ambulance, EMS, Fire
- Long Term Care

### FEATURES

- 18 or 15 UVC emitters
- External UVC sensors
- LED operations indicator
- Emergency stop switch
- Door hanger warning signs
- Treatment tracking tablet
- Dual-tech motion sensors
- Remote control
- Wide, anti-tip casters

### SPECIFICATIONS

- Size: 60.5"H x 28.5"W x 27.5"D
- Weight: 126 lbs.
- Construction: Aluminum/  
Stainless Steel
- Electrical: 20A @ 120V / 10A @  
240V

99.99% effective at killing  
viruses and bacteria



**CREST**  
HEALTHCARE®

800.328.8908

800.369.9207

CrestHealthcare.com