

UVC Safezone

Maximize Power. Maximize Mobility.

Adopting UVC Cleaning Systems disinfection technology will complement your existing workflow processes without requiring additional full-time employees. Our products allow your current staff to **work smarter, not harder** with easy to operate, tested and proven, medical grade, portable UVC Disinfection devices.

UV EFFECT of 254nm WAVELENGTH

Photochemical destruction of the microorganisms' DNA by absorption of UV-C light which shuts down the organism's reproductive processes rendering it non-pathogenic.

FEATURES

- Wireless on/off remote
- (9) 95W UV lamps
- 30 min adjustable UV timer
- Run 2 units with 15 Amp plug

APPLICATIONS

- Hospitals, Clinics, Medical
- Schools, Commercial
- Home, Residential
- Ambulance, EMS, Fire
- Long Term Care

SPECIFICATIONS

- **Size without Base:** 31"H x 13"W x 13"D
- **Size with Base:** 45"H x 18"W x 18"D
- **Construction:** Aluminum/Stainless Steel
- **Electrical:** 7.5A @ 120V / 4A @ 240V



*99.99% effective at killing
viruses and bacteria*

CREST
HEALTHCARE®



800.328.8908



800.369.9207



CrestHealthcare.com