

self-locking rocking chairs

# Thera-Glide®

R, W and T Series

FOCUSED ON WHAT MATTERS



SPECIALTY  
HEALTHCARE  
SEATING



WOOD MODEL (W)



METAL MODEL (R)



TRANSITIONAL  
MODEL (T)

MAINTAINING PATIENT AUTONOMY

*Optima*  
PRODUCTS INC.

# Thera-Glide®

Self-locking rocking chairs


## R Series Decorative Metal

(Widths available: 17", 20" and 24")

The Thera-Glide® glider has a place for itself in the healthcare world, both for its therapeutic effects and for its incredible functional and safety features. Its patented automatic immobilization mechanism keeps the chair fixed when getting up and allows rocking when the person is fully seated. The system is easy, quiet and locks in any position. Depending on the model, its adjustment options and accessories, you can adapt the Thera-Glide® to meet your specific needs. It is easy to maintain thanks to the removable upholstery and the materials used in its design.



Secure when sitting down



Rocks freely



Secure when getting up



Contemporary style backrest for easy cleaning (Classic and Comfort Plus styles also available)

Elongated armrests for a better grip

Seat depth and backrest angle adjustment standard on 17" and 24" models, optional on the 20" Metal model

Fast and easy indexer height adjustments

With a metal frame with a lifetime warranty, a contemporary style backrest with smooth surfaces and the polyurethane armrest option, the Thera-Glide model is a safe rocking chair that meets the strictest disinfection standards



# W Series Wood

(Widths available: 17", 20" and 24")



CLASSIC BACKREST



COMFORT+ BACKREST



Height-adjustment  
is standard for  
optimal comfort

Seat depth and backrest angle adjustment  
standard on 17" and 24" models,  
optional on the 20" Wood model

The Thera-Glide®'s wood components  
are attached to a robust and durable  
inner steel frame.

5 wood stains and three backrest  
styles enable you to customize  
the Wood Thera-Glide® to your liking.



10. CHERRY  
(R/W/T Series)



11. NATURAL  
(W/T Series)



12. MAHOGANY  
(W/T Series)



13. WALNUT  
(W/T Series)



14. EBONY  
(W/T Series)



Lateral stabilizers for increased safety

# T Series Transitional Style

(Single width: 22")

The Thera-Glide® Transitional provides all the safety features of our other rocking chairs. Its "transitional" design blends comfort, warmth and modernity and will prove to be a valuable asset for your interiors.

Our T Series is suitable for many people. Thanks to its easy-to-clean waterfall cushion, curved wooden armrests, adjustable height and maximum capacity of 350 lbs. (159 kg), the Thera-Glide® Transitional is a safe rocking chair ideal for the common areas of your residences.



Ornamental wood frame  
(see choice of W Series colours)

Upholstered side panels  
for a sleek, modern finish

Rocking mechanism manual  
lock (standard on all models)





# Rocking...

- + Significantly improves the quality of life for people experiencing a loss of autonomy.
- + Leads to a reduction in medication intake.
- + Therapeutically combats anxiety and depression in people with dementia and/or Alzheimer's<sup>1</sup>.
- + Provides significant benefits to newborns and their mothers<sup>2</sup>.

<sup>1</sup>[Child Study Center, School of Medicine, Yale University]. <sup>2</sup>[Watson, Wells and Cox, Rochester University, 1998].

## BASIC FEATURES

Thera-Glide® chairs come standard with vinyl upholstery. Silververtex vinyl containing silver nitrate\*, Cryptons as well as an almost infinite selection of coverings are available as options. Our customers can also provide us with their own fabric if needed. (COM).

\*Silver nitrate is known for its antibacterial properties.

## SERIES SPECIFICATIONS

DESCRIPTION	R Series	W Series	T Series
Extended armrests with ergonomic palm grips	+	+	+
Removable upholstery for easy cleaning	+	+	+
Automatic infinite position lock	+	+	+
Glide motion lock-out	+	+	+
Anti-tip side stabilizers	+	+	+
Rear wheels to aid moving	+	0	0
Adjustable seat height	+	+	+
Sinuous loop spring seat constructions	—	—	+
Matching colour undercarriage cover	—	+	+
+: Standard    0: Optional    —: Not available			

## OPTIONS AND ACCESSORIES

ITEM	DESCRIPTION	R Series	W Series	T Series
R-C1	Adjustable seat depth and backrest angle 15" to 18" (38 to 46 cm) / 95 to 115 degrees	0	0*	—
R-C2	Upholstered side panels	0	—	—
R-C4	Adjustable height kit for 21" to 24"	0	—	—
R-C7	Polyurethane armrests	0	—	—
R-020	Gliding metal ottoman	0	—	—
R-030	Fixed metal ottoman	0	—	—
RA-01	Opaque activity tray	0	0	—
RA-02	Activity tray with Lexan window	0	0	—
RA-03	Lateral tray straps	0	0	—
RA-03A	Wrap-around tray straps	0	0	—
RA-04	Rear wheels	—	0	0
RA-09	Side storage bag	0	0	—
RA-10	1" (2.5 cm) seat insert	0	0	—
RA-13	Seat wedge insert	0	0	—
RA-19	Universal perineal abductor kit with black fixing mat	0	0	—
RA-23	4" (10 cm) lateral trunk supports with backrest canvas	0	0	—
RA-24	6" (15 cm) lateral trunk supports with backrest canvas	0	0	—
RA-28-D	Undercarriage cover delete	—	0	—
RA-29	Backrest switch delete	—	0	0
RB-18	Cervical support	0	0	—
WOF-22	Upholstered wood ottoman with footboard	—	0	0
Z-C8	PRE-CURE® PRO seat cushion	0	0	—
+: Standard    0: Optional    —: Not available				

\* Optional on 20" and standard on 17" and 24".

## TECHNICAL SPECIFICATIONS

	R Series – Decorative Metal			W Series - Wood			T Series - Transitional Style
	R-505 SMALL	R-605 REGULAR	R-705 LARGE	W-2XXX-17-XX SMALL	W-2XXX-20-XX REGULAR	W-2XXX-24-XX LARGE	
Total height	44" to 47" / 112 cm to 119 cm			44" to 47" / 112 to 119 cm			42 ½" to 45 ½" / 108 to 116 cm
Seat-to-floor height	Adjustable from 16" to 19" / 41 to 48 cm						
Total width	22" / 56 cm	26" / 66 cm	29" / 74 cm	22 ½" / 57 cm	26 ½" / 67 cm	29 ½" / 75 cm	27 ¼" / 69 cm
Seat width	16 ½" / 42 cm	20 ½" / 52 cm	23 ½" / 60 cm	16 ½" / 42 cm	20 ½" / 52 cm	23 ½" / 60 cm	22" / 56 cm
Seat-to-armrest height	8" / 20 cm			8" / 20 cm			7" / 18 cm
Depth	32" / 81 cm			33" / 84 cm			33" / 84 cm
Seat depth	17 ½" / 45 cm			17" / 43 cm			18" / 46 cm
Backrest height	31" / 79 cm			29" / 74 cm			28" / 71 cm
Maximum weight capacity	300 lbs. / 136 kg		450 lbs. / 205 kg	275 lbs. / 125 kg		450 lbs. / 205 kg	350 lbs. / 159 kg
Weight	58 lbs. / 26 kg	61 lbs. / 28 kg	67 lbs. / 31 kg	57 lbs. / 26 kg	60 lbs. / 27 kg	66 lbs. / 30 kg	75 lbs. / 34 kg

Please note that the technical specifications are subject to change without notice.

# our complete product range

## MAINTAINING PATIENT AUTONOMY

THERA-GLIDE®  
WOOD



THERA-GLIDE®  
METAL



THERA-GLIDE®  
TRANSITIONAL



DYN-ERGO



NOMAD



VERSATILE



EVOLUTION



UTILITY



## COMFORT AND POSITIONING

COMP-POSITION



HARMONY



OPTIMUM



## ACCESSORIES

PRE-CURE® PRO  
CUSHION



ANTI-SLIP  
SURFACES



Visit our website or contact us for more details about our range of products.

lpa+

[lpamedical.com](http://lpamedical.com)

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA  
T 418-681-1313 1-800-663-4863 F 418-681-4488

*Optima*  
PRODUCTS INC.

2230 Edgewood Ave S. Suite 3A Minneapolis, MN 55426  
952-745-0027 | 888-990-9394 | [optimaol.com](http://optimaol.com)