



Where seating is a science.



OPTIMA LIFT CHAIRS & RECLINERS



L-1012-ZS Wall-Saver Lift Chair Transitional
Style w/ Plush Back
Healthcare Vinyl



R-1001-ZS Wall-Saver Recliner
2 Pillow Back w/ EZ Reach Handle
Crypton & Naugasoft

DESIGNED SPECIFICALLY FOR HEALTH CARE FACILITIES

All the features you need in an institutional grade lift chair and recliner including contract grade moisture impermeable fabrics, wall-saver designs, removable seat cushion for easy cleaning, replaceable parts, and more.



OPTIMA HEALTH CARE LIFT CHAIRS

Everything about the Optima Lift Chair is made with today's busy healthcare facility in mind. Its removable chaise cushion and moisture-impermeable fabric make the Optima Lift Chair easy to clean and disinfect. Most models are available in a wall-saver design to prevent wall damage and fit where ordinary lift chairs won't.

Available in 3 styles: Transitional, Transitional w/Plush Back, 2 Pillow Back
Available Sizes: Petite, Standard, Wide, & X-Wide



The hand-held control permits left hand, right hand, or attendant operations. A variety of contract grade vinyl and Crypton fabrics complement any décor. Customer's Own Material (COM) welcome. Safety and durability come first with a full chaise leg-rest covering the recliner mechanism for a pinch-free design, a wide extra stable all steel base, heavy duty mechanisms, hefty lifting capacities, TB-117 and NFPA compliant construction, replaceable parts, and a lifetime limited warranty on the lift mechanism. Optima Lift Chairs meet your facility's specifications. Wall-Saver design requires only 4 inches of wall clearance, uses less floor space and protects walls from costly, unsightly damage.



LIFT CHAIR OPTIONS:

Heat/Massage, Holster for Hand Control, Side Pocket Right, Left, or Both, Cervical Pillow, Fabric Options
Available Sizes - Petite, Standard, Wide, X-Wide

OPTIMA HEALTH CARE RECLINER

Cleanability - Impermeable fabric covers all surfaces, even underneath the removable chaise cushion.

Safety - Chaise leg rest covers recliner mechanism. No pinch points, wide stable base.

Wall-Saver Design - Only 4 inches of clearance needed. Reduces damage to walls & chair backs while conserving space.

Crypton, Leatherette, or COM fabrics - We'll source, grade, and provide designer specific fabrics.

Exceptional Durability - Heavy duty recliner mechanism and hardwood construction.

Fire Codes Compliance - Meets CA 117 and NFPA Life Safety Codes.

Range of Styles & Sizes - Available in Transitional, Transitional w/ Plush Back, Petite, Standard, Wide and X Wide.



R-1001-ZS Wall-Saver Manual Recliner
2 Pillow Back w/EZ Reach Handle



R-1012-ZS Wall-Saver Manual Recliner Transitional Style
w/ Plush Back & Standard Handle

Manual Recliner available in 3 different Ottoman Spring Tensions: Light, Medium, Firm



R-1002-ZS Wall-Saver Manual Recliner
Transitional Style w/ Standard Handle

**Power Recline also available:
PR-1001, PR-1002, PR-1012**

Complete specifications available upon request or online at www.optimaol.com/lfrclnspecifications

RECLINER OPTIONS:

Standard Handle, EZ Reach Handle, Ottoman Spring Tension, Power Recline, Heat/Massage (Power Only), Side Pocket Right, Left, or Both, Cervical Pillow, Fabric Options
Available Sizes - Petite, Standard & Wide

OPTIONS

Heat & Massage, Side Pockets, Cervical Support, Recliner Handle, Stock Fabric



Heat Massage Hand Control
several different massage
configurations available



Hand Control Holster



Side Pocket
available for right,
left or both sides



Cervical Head Rest



Standard Recliner Handle



EZ Reach Handle

STOCK FABRICS



Naugasoft Cocoa




Naugasoft
Cabernet/Burgundy



Naugasoft Harness



Naugasoft Sapphire Blue



Naugasoft Deep Sea Green

Crypton® Stain Resistant fabrics and Reduced Environmental impact fabrics available.

Complete specifications available upon request or online at www.optimaol.com/lfrclnspecifications

OUR PRODUCTS



Transitional Style
Thera-Glide
Safety Glider



Wood
Thera-Glide
Contemporary
Style



Wood
Thera-Glide
Traditional
Mission Style



Dyn-Ergo
Scoot Chair



Evolution
Mobility Chair



Optimum
Posturo-Pedic
Positioning Chair



Comp-Position
Positioning Chair



www.optimaol.com

952-745-0027

Optima Products Inc.

2230 Edgewood Avenue | Suite 3A
St. Louis Park, MN 55426

Toll free: 888-990-9394 | Fax: 952-745-0197
Email: sales@optimaol.com

