

SAFE, STABLE AND COMFORTABLE MOBILITY



PROTEKT® PILOT UPRIGHT WALKER

FEATURES:

- Adjustable Height Armrests**
 Deluxe padded armrests enable users of a wide range of heights to stand erect. This feature allows users to look ahead, so they feel safe, secure, and confident while walking.
- Lockable Handbrakes**
 Ergonomic dual handbrakes provide better control for braking, sitting, and parking.
- Comfortable Seat**
 Deluxe fabric seat provides a sturdy and safe seating surface.
- Backrest Support**
 Detachable wide backrest provides comfortable and secure back support.
- Stability & Security**
 8" front wheels swivel 360° for easy and safe maneuvering. Strong enough for indoor, outdoor, and uneven surfaces.
- Stand Assist Handles**
 Comfortable padded handles provide an easier and safer transition between sitting and standing.
- Compact Design**
 Allows for easy transport and convenient storage.
- Convenience Package Included**
 Beverage holder, safety LED light, personal item bag and cane holder.

ISSUES:



OUR SOLUTION:



1. Slouching or hunching over.	1. Allows users to stand up straight with improved posture.
2. Unbalanced/anxiety of falling.	2. Users feel safer and more secure with added stability.
3. Uncomfortable and causes pain.	3. More comfortable positioning mitigates pain in the back, hands, and wrists.
4. Reduction of independence and freedom.	4. Allows users to be more mobile, active, and walk longer distances with less stress.

PROTEKT® PILOT

is designed with simple features for ease of use.

- | | |
|---------------------------|-----------------------|
| 1. Handgrip | 9. Front wheel |
| 2. Parking/handbrake | 10. Seat |
| 3. Armrest pad | 11. Seat rail |
| 4. Height adjustment tube | 12. Backrest |
| 5. Height adjustment lock | 13. Personal item bag |
| 6. Stand assist handles | 14. Beverage holder |
| 7. Rear wheel | 15. Safety LED light |
| 8. Brake adjustment knob | 16. Cane holder |



STAND ASSIST HANDLES



FOLDS EASILY FOR TRANSPORT



COMPACT DESIGN FOR STORAGE

Features	Specifications
User Height Range	5'0" - 5'11"
User Weight Capacity	300 lbs.
Dimensions	29.6"Lx25"Wx41"H
Armrest Height Range	41" - 46"
Seat Dimensions	19"Wx10"D
Seat Height from Floor	21"
Rear Wheel Base	25"
Wheel Size	8" diameter
Unit Weight	22 lbs.
Folded Dimensions	29.6"Lx9.5"W
Warranty	1 Year Limited

Ordering Information	
UPWLKALDLX19BL	Protekt® Pilot - BLUE
UPWLKALDLX19RD	Protekt® Pilot - RED
UPWLKALDLX19SL	Protekt® Pilot - SILVER

For more information, please call or visit:

270 Washington Street,
Mount Vernon, NY 10553

Tel: 855-237-7622

www.proactivemedical.com