

# Prevent Falls & Wandering

## CordLess® Pads & Mats For Bed & Chair Exit Alarms



Works with:



Central  
Monitor Unit  
433-CMU-40/60



CordLess®  
Alarm  
433-EC



CordLess®  
Monitor  
TL-2100G

### No Cords Between Pads & Monitors

When using a CordLess® Pad with a CordLess® Monitor you are able to remove the alarm noise from the room, creating a quieter environment for residents.

The Pad is placed under resident, either on a chair or bed, and the floor mat can be placed either by the bed or by the door. When a resident gets up from a pad, or steps on a floor mat, a wireless signal is sent to the monitor setting off an alarm to alert caregivers.

MFR#	SMART PRODUCT DESCRIPTION
GBT-WI	CordLess® Bed Pad 20" x 30" (1 year pad life)
GBT-RI	CordLess® Bed Pad 10" x 30" (1 year pad life)
GBT-90, GBT-45	CordLess® Bed Pad 10" x 30" (90 day, 45 day pad life)
GCT-WI	CordLess® Chair Pad 10" x 15" (1 year pad life)
GCT-90, GCT-45	CordLess® Chair Pad 10" x 15" (90 day, 45 day pad life)
FMT-03C	CordLess® Floor Mat 24" x 36" Black (1 year warranty)
FMT-05C	CordLess® Floor Mat 24" x 36" Gray (1 year warranty)
FMT-07C	CordLess® Floor Mat 24" x 48" Gray (1 year warranty)
GWC-02	CordLess® Fall Mattress 70.5" x 29.5" x 1.5" (90 day warranty)
LM-01C	CordLess® Landing Mat 24" x 71" x 1" Gray (1 year warranty)



Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)