

DeVilbiss 5 Liter Oxygen Concentrator

HCPCS: E1390



FIGURE A



FIGURE B

- Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter Oxygen Concentrator was designed with enhanced quality, durability and simplicity in mind
- Patented DeVilbiss Turn-Down Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator
- Exclusive DeVilbiss OSD (Oxygen Sensing Device), standard on every unit, ensures patient safety and reliability for longer service intervals
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device
- Sold separately, the Transfill Caddy is designed with a steel frame to simply connect the 18 x 14 inch shelf, designed to support up to 40 lbs, to the DeVilbiss 5L Concentrator. Rear view (Figure A), Front View (Figure B)

PRODUCTS

Description

525DS	DeVilbiss 5 Liter Oxygen Concentrator, EA
525DS-EW	DeVilbiss 5 Liter Oxygen Concentrator, EA
525DS-Q	DeVilbiss 5 Liter Oxygen Concentrator, EA

SPECIFICATIONS

Description

Delivery Rate	0.5 to 5 LPM
Outlet Pressure	8.5 psig
Oxygen Percentage	1-5 LPM, 87% - 96%
Power Consumption	310 Watt Average, 275 Watts @ 1.2 LPM and Below
Product Dimensions	24.5" (H) x 13.5" (W) x 12" (D)
Product Weight	36 lbs.
Warranty	3 Years

ACCESSORIES

Description

515LF-607	FLOW METER PACKAGE, LOW OUTPUT, EA
525DD-650	Transfill Rack System, 525DS Series, 1/e, EA
MQ0412	No Smoking Sign, 25/CV, CS
PF1100TUB	Transfill Connection Tubing, 1/ea, EA

PARTS

Description

525DD-606	Concentrator (5L, 10L) Outlet Port, 1/ea, EA
525DD-603	Base Assembly, 1/ea, EA
444-566	LADDER CLAMP - DFT, EA
PV5LD-651	FINAL BACTERIA FILTER - DFT, EA
505DZ-645	POWER CORD STRAIN RELIEF, EA

PARTS	Description
525DD-639	DOOR INTAKE FILTER PACKAGE, EA
525DD-602	REAR COVER PACKAGE, EA
303DZ-628	CABINET SCREWS 6-PKG, EA
525DD-628	CABINET SCREW PACKAGE, EA
525D-610	ACC. TANK S.S.P., EA
515ADZ-627	FUSE .125 A FOR 515ADZ/S-425, EA
PVO2D-607	PURGE HARN/CK VALVE 2-PK, PKG
515ADZ-702	ROTARY VALVE REPLACEMENT, EA
525D-608	COMMUNICATION WIRE HARNESS, EA
515DZ-634	COOLING FAN PACKAGE, EA
515ADZ-633	FAN GUARD, EA
525DD-614	INTAKE CANISTER PACKAGE, EA
515A-705	EXHAUST MUFFLER PKG, EA
525DD-605	COMPRESSOR BOX PACKAGE, EA
505DZ-634	SILICON TUBING (2'), EA